

Manufacturer and Supplier of Laboratory Products

V.2022

Biologix Group Ltd

Biologix is specialized in the R & D and production of disposable medical detection reagents & supplies, laboratory supplies & instruments, and complete biobanking solutions. Biologix has a responsible and flexible team of major technological and managerial talents in the field of life science and medicine. Biologix is rated as a High-Tech Enterprise. Taking "Serve Human Health" as the object, Biologix's two brands: BIOLOGIX and CryoKING enjoys international excellence and reputation from the power of the sustainable innovation and continuous improvement.

Biologix focuses on the R & D and innovation of core products, and has the standard class 100,000 clean room which are compliant with GMP standards and number of professional R & D testing laboratories which are all ISO 9001, CE, and FDA registered, and follow the strict CDC guidelines for purity and efficiency. Biologix is committed to the design and manufacture of high-quality products, enriching product lines, and strengthening the service of complete solutions.

Biologix has developed a focus on constant innovation and improvement to offer optimal products, innovative solutions, and service to about dozens of countries and regions in Europe, America, and Asia. Biologix has been an international and well-known company in the world. Biologix wishes to strengthen the cooperation with international companies in OEM service and technical cooperation hand in hand to contribute to human health.

Ordering and Services



Biologix offers a Private Labeling Program, providing high quality, branded lab supplies at very competitive prices. Many Biologix products are currently private labeled and sold worldwide by some of the most prominent brands and distributors.



Biologix also operates as an OEM provider for plastic laboratory consumables. With through our inhouse molding and tooling capabilities, we can provide the exact products needed by our customers at significantly lower costs than other companies can achieve. Biologix currently supplies some of the world's top distributors with its OEM products.



A. Free Product Samples

Let us send you a FREE sample package of our great products. Contact us today if you are interested in trying out Biologix products before you buy.

B. Distribution

If interested in distribution opportunities, please contact us for distributor terms and pricing. We are accepting distribution opportunities worldwide. Please contact us at *Inquiry@BiologixGroup.com*.

Guaranteed Quality Exceptional Service

Professional Design Logical Pricing











CE ISO FDA

Contents

| General Laboratory Sup | pplies | Cell Culture Plates | |
|-------------------------------|--------|--------------------------------|-------|
| - | | Serological Pipettes———— | |
| Centrifuge Tubes———— | | Cell Culture Slides———— | |
| Microcentrifuge Tubes———— | | Cell Strainers————— | |
| Centrifuge Tube Racks———— | | Sterile Syringe Filters & Vacu | |
| Pipet Tips—————— | | Cell Scrapers & Lifters——— | |
| Filter Tips—————— | | | |
| Rainin LTS Pipette Tips | | PCR | |
| Automation Tips————— | | PCR Tubes | 50 |
| Transfer Pipets | | 8-Strip PCR Tubes———— | |
| Solution Basins————— | | PCR Plates————— | |
| Beakers | | Terriaces | 01,02 |
| Reagent Bottles | | | |
| Test Tubes—————— | | Microbiology | |
| Slide Storage Boxes————— | | Microbiology | |
| Deep Well Plates | | Petri Dishes | 64 |
| Sealing Films————— | | Inoculating Loops & Needles | 65 |
| Gloves | 29 | Cell Spreaders———— | |
| CryoKING Biobanking | I | | |
| Screw Cap Microtubes | 31 | Drosophila | |
| Cryogenic Vials | 32 | • | |
| Traditional Cryogenic Vials | 33 | Drosophila Vials———— | |
| CryoKING Cryogenic Vials———- | 34 | Drosophila Blugs | |
| CryoKING Combo | | Drosophila Plugs | /0 |
| Automated Product Series of | | | |
| SBS Standards————— | 36 37 | | |
| Premium Cardboard Cryogenic B | | Sample Collection | and |
| | | Detection Products | |
| ID-Color™ Cardboard Cryogenic | | Sample Collection | |
| | | and Detection Products——— | 71 |
| PC Cryogenic Boxes | 40 | and Detection Froducts | 71 |
| PP Cryogenic Boxes | 41 | | |
| Freezer Racks | 42,43 | | |
| Screw Cap Capper/Decapper&Sca | | Instruments | |
| Tube Handling System | • • | Liquid Handling | 74,75 |
| Biobanking Management Softwa | | Collection System | |
| low Temperature Freezers———— | | Stirrers————— | |
| Liquid Nitrogen Tanks————— | | Mixers & Rotators | |
| Liquid Millogell Idliks | 40 | Shakers | 79 |
| | | Thermocyclers————— | 80 |
| Cell Culture Series | | Centrifuges | |
| | | Dry Bath Incubators———— | |
| Cell Culture Flasks————— | | Electrophoresis————— | |
| Cell Culture Dishes | 51 | BluView Transilluminator—— | 88,89 |

Catalog Number Index

| 01-110X 78 02-8024 27 07-3060 51 20-1250 10 01-120X 78 02-6022 27 07-3150 51 21-0010 10 01-2101 74 02-8005 27 07-3001 53 21-0200 10 01-2101 74 02-8005 27 07-5001 53 21-0200 10 01-2104 74 03-1300 87 07-5005 53 21-1280 10 01-2115 74 03-2300 87 07-5005 53 21-1280 10 01-2116 74 03-300 87 07-5005 53 21-1280 10 01-2115 74 03-3210 87 07-5055 53 23-0011 11 01-2115 74 03-3307 87 07-5055 53 23-0011 11 01-2112 74 03-3400 87 07-5056 53 23-0011 11 01-2122 74 03-4221 85 07-6006 52 23-0020 11 01-2120X 75 03-6210 88 07-6024 52 23-0021 11 01-220X 75 03-6210 88 07-6024 52 23-0001 11 01-221X 74 04-0050 24 07-6048 52 23-0001 11 01-221X 74 04-0060 24 07-6096 52 23-0011 51 01-221X 74 04-0060 24 07-8055 50 23-00515 12 01-241X 74 04-01250 24 07-8057 50 23-00515 12 01-243X 74 04-01250 24 07-8057 50 23-00515 12 01-243X 74 04-01250 24 07-8057 50 23-00515 12 01-243X 74 04-01250 24 07-8057 50 23-00515 12 01-240X 77 04-10500 24 07-9075 50 23-00515 12 01-2801 76 04-05000 24 07-9075 50 23-00515 12 01-2803 74 04-10600 24 07-9075 50 23-00515 12 01-2801 76 04-05000 24 07-9075 50 23-00515 12 01-2803 74 04-10500 24 07-9075 50 23-00515 12 01-2803 77 04-10500 24 08-PES322R 56 25-0052 20 01-330X 77 04-10500 24 08-PES322R 56 25-0052 20 01-330X 77 04-10500 24 08-PES322R 56 25-0052 20 01-330X 77 04-10500 24 08-PES322R 56 25-1005 20 01-330X 77 04-10500 24 08-PES322R 56 25-1005 20 01-330X 77 04-10500 24 08-PES3250 56 30-013841 21 01-20X 89 04-5000 24 08-PES3250 56 30-013841 21 01-90X 81 04-2060 24 08-PES3250 56 30-013841 21 01-90X 81 04-2060 24 08-PES3250 56 30-013841 21 01-90X 81 04-2060 24 08-PES3250 56 30-013841 21 01-90X 80 04-2050 24 10-9151 07 51-0501 68 01-90X 81 04-2060 24 10-9151 07 51-0501 68 01-90X 81 04-2060 24 10-9151 07 51-0501 68 01-90X 80 04-2050 24 10-9151 07 51-0501 68 01-90X 81 04-2060 24 10-9151 07 51-0501 68 01-90X 81 04-2060 24 10-9151 07 51-0501 68 01-90X | Cat # | Page# | Cat # | Page# | Cat # | Page# | Cat # | Page# |
|--|----------|-------|---------|-------|-------------|-------|-----------|-------|
| 01-190X 78 02-7022 27 07-3150 51 21-0010 10 01-2102 74 02-9010 27 07-5002 53 21-1000 10 01-2104 74 03-300 87 07-5002 53 21-1000 10 01-2104 74 03-300 87 07-5002 53 21-1000 10 01-2115 74 03-3070 87 07-5025 53 23-0010 11 01-2115 74 03-3070 87 07-5025 53 23-0011 11 01-2115 74 03-3070 87 07-5025 53 23-0011 11 01-2122 74 03-4221 85 07-6006 52 23-0001 11 01-2122 74 03-4221 85 07-6006 52 23-0001 11 01-2212 74 03-4221 85 07-6006 52 23-0001 11 01-220X 75 03-6210 88 07-6024 52 23-1000 11 01-220X 75 03-6210 88 07-6024 52 23-1001 11 01-221X 76 04-0008U 24 07-6096 52 23-0011 12 01-221X 74 04-0030U 24 07-6096 52 23-0011 12 01-241X 74 04-0030U 24 07-6096 52 23-00115 12 01-241X 74 04-0030U 24 07-8026 50 23-00115 12 01-244X 74 04-015U 24 07-8076 50 23-00115 12 01-244X 74 04-015U 24 07-8076 50 23-00115 12 01-244X 74 04-0250U 24 07-8076 50 23-00115 12 01-220X 76 04-0500U 24 07-8076 50 23-00115 12 01-220X 76 04-0500U 24 07-8076 50 23-00115 12 01-243X 74 04-0030U 24 07-8076 50 23-00115 12 01-243X 74 04-0050U 24 07-8076 50 23-00115 12 01-243X 74 04-0050U 24 07-8076 50 23-00115 12 01-2801 76 04-0500U 24 07-8075 50 23-00115 12 01-2803 74 04-1090U 24 08-PES332S 56 25-0051 20 01-330X 77 04-1050U 24 08-PES332S 56 25-0051 20 01-330X 77 04-1050U 24 08-PES332S 56 25-1025 20 01-330X 77 04-1250U 24 08-PES32S 56 25-1025 20 01-330X 77 04-1250U 24 08-PES32S 56 25-1025 20 01-330X 80 04-2015 24 08-PES32S 56 25-1025 20 01-330X 80 04-2050 24 08-2050 27 51-1050 26 | 01-110X | 78 | 02-3024 | 27 | 07-3060 | 51 | 20-1250 | 10 |
| 01-2101 74 02-8005 27 07-5001 53 21-0200 10 01-2102 74 03-3100 87 07-5005 53 21-1200 10 01-2104 74 03-1300 87 07-5005 53 21-1200 10 01-2110 74 03-2310 87 07-5010 53 23-0010 11 01-2112 74 03-3070 87 07-5055 53 23-0010 11 01-2112 74 03-3100 87 07-5055 53 23-0010 11 01-2115 74 03-322 85 07-6012 52 23-1000 11 01-2155 74 03-422 85 07-6012 52 23-1000 11 01-2155 74 03-422 85 07-6012 52 23-1000 11 01-221X 76 03-6021 88 07-6024 52 23-0010 11 01-221X 76 04-0008U 24 07-6048 52 23-0010S 12 01-241X 74 04-0030U 24 07-6048 52 23-0010S 12 01-243X 74 04-0030U 24 07-6056 52 23-0010S 12 01-243X 74 04-0050U 24 07-8075 50 23-00215 12 01-2270X 76 04-0500U 24 07-9075 50 23-00205 12 01-2801 76 04-0500U 24 07-9075 50 23-00005 12 01-2803 74 04-1010U 24 07-9075 50 23-00005 12 01-2803 74 04-1010U 24 07-9175 50 23-0101S 12 01-2803 74 04-1016U 24 07-9175 50 23-0101S 12 01-2803 74 04-1016U 24 07-9175 50 23-0101S 12 01-2803 74 04-1016U 24 08-PES322R 56 25-0025 20 01-320X 77 04-1030U 24 08-PES322R 56 25-0052 20 01-330X 77 04-1050U 24 08-PES3322R 56 25-0052 20 01-330X 80 04-2015 24 08-PES3322R 56 25-0052 20 01-330X 80 04-2015 24 08-PES3322R 56 25-0052 20 01-330X 80 04-2055 24 09-PES2250 56 30-0338A1 21 01-700X 80 04-2056 24 09-PES2250 56 30-0338A1 21 01-700X 80 04-2056 24 09-PES2250 56 30-0388A1 21 01-700X 80 04-2056 24 10-9515 07 51-0600 68 01-960X 86 04-3500 24 10-9515 07 51-0600 68 01-960X 86 04-3500 24 10-9515 07 51-0600 69 01-960X 86 04-2550 24 10-9515 07 51-0600 69 01-96 | 01-120X | 78 | 02-6022 | 27 | 07-3100 | 51 | 20-1251 | 10 |
| 01-2102 74 02-9010 27 07-5002 53 21-1000 10 01-2104 74 03-3030 87 07-5005 53 21-1000 10 01-2115 74 03-2310 87 07-5010 53 23-0010 11 01-2115 74 03-3070 87 07-5025 53 23-0011 11 01-2115 74 03-3100 87 07-5025 53 23-0011 11 01-2122 74 03-4221 85 07-6006 52 23-0201 11 01-2122 74 03-4221 85 07-6012 52 23-1001 11 01-220X 75 03-6210 88 07-6024 52 23-1001 11 01-221X 76 04-00088 24 07-6024 52 23-1001 11 01-221X 76 04-0015U 24 07-6096 52 23-0011 11 01-221X 76 04-0015U 24 07-6096 52 23-0011 12 01-242X 74 04-0030U 24 07-8075 50 23-0011S 12 01-242X 74 04-0050U 24 07-8075 50 23-0011S 12 01-244X 74 04-015U 24 07-8075 50 23-0011S 12 01-244X 74 04-015U 24 07-8075 50 23-0011S 12 01-244X 74 04-0250U 24 07-8175 50 23-0101S 12 01-2405 76 04-0600U 24 07-9075 50 23-0201S 12 01-241X 74 04-0100 24 07-9075 50 23-0201S 12 01-2805 74 04-1008U 24 07-9175 50 23-0101S 12 01-2805 74 04-1015U 24 07-9175 50 23-0101S 12 01-2805 74 04-1015U 24 08-PES1322R 66 25-0051 20 01-330X 77 04-1030U 24 08-PES322R 66 25-0051 20 01-330X 77 04-1030U 24 08-PES322R 66 25-0051 20 01-330X 77 04-1030U 24 08-PES3328 66 25-1025 20 01-330X 77 04-1030U 24 08-PES3328 66 25-1025 20 01-330X 77 04-1030U 24 08-PES3328 66 30-0135 21 01-261XA 79 04-1010U 24 08-PES3328 66 30-0135 21 01-270X 80 04-2015 24 08-PES3328 66 30-0135 21 01-281X 79 04-1000U 24 08-PES3328 66 30-0135 21 01-300X 80 04-201S 24 08-PES3328 66 30-0138 21 01-700X 80 04-2015 24 08-PES3250 56 30-0138 21 01-700X 80 04-2015 24 09-CA2250F 56 41-6005 26 01-830X 82 04-215 24 09-CA250F 56 41-6005 26 01-830X 83 04-2500 24 10-9510 07 51-0501 68 01-900X 86 04-3015 24 10-9510 07 51-0501 68 0 | 01-130X | 78 | 02-7022 | 27 | 07-3150 | 51 | 21-0010 | 10 |
| 01-2014 74 03-1300 87 07-5005 53 21-12550 10 11-2110 74 03-2310 87 07-5010 53 23-0010 11 11 01-2115 74 03-3070 87 07-5055 53 23-0010 11 11 01-2112 74 03-3100 87 07-5055 53 23-0011 11 11 01-2112 74 03-3100 87 07-5056 53 23-00201 11 11 01-2155 74 03-4221 85 07-6006 52 23-0201 11 01-2155 74 03-4222 85 07-6012 52 23-1000 11 10 11 01-2121X 76 03-6210 88 07-6024 52 23-1001 11 01-2121X 76 04-0008U 24 07-6048 52 23-0010S 12 01-221X 77 04-00060 24 07-6048 52 23-0010S 12 01-241X 74 04-0030U 24 07-6048 52 23-0010S 12 01-241X 74 04-0030U 24 07-8075 50 23-0011S 12 01-243X 74 04-0060U 24 07-8075 50 23-0021S 12 01-243X 74 04-0050U 24 07-8075 50 23-0051S 12 01-243X 74 04-0050U 24 07-8075 50 23-0021S 12 01-243X 74 04-0050U 24 07-8075 50 23-0021S 12 01-2805 74 04-0050U 24 07-8075 50 23-0021S 12 01-2805 74 04-1010U 24 07-9075 50 23-0020S 12 01-2805 74 04-1010U 24 07-9075 50 23-0020S 12 01-2805 74 04-1010U 24 08-PES1322R 56 25-0051 20 01-2803 74 04-1015U 24 08-PES3222R 56 25-0051 20 01-330X 77 04-1060U 24 08-PES3322R 56 25-0052 20 01-330X 77 04-1060U 24 08-PES3322R 56 25-1025 20 01-330X 77 04-1060U 24 08-PES332R 56 30-0138 21 01-700X 80 04-2028 24 08-PES332P 56 30-0138 21 01-700X 80 04-2028 24 08-PES332P 56 30-0138 21 01-700X 80 04-2028 24 08-PES332P 56 30-0138 21 01-700X 80 04-2030 24 09-PES2505 56 41-51025 26 01-830X 82 04-2105 24 09-02250 56 41-5100 26 01-830X 82 04-2105 24 09-02250 56 41-5100 26 01-830X 82 04-2105 24 09-02250 56 41-5100 26 01-830X 82 04-2250 24 10-0152 07 51-0500 68 01-900X 86 04-3008 24 10-9050 07 51-0500 68 01-900X 86 04-3008 2 | 01-2101 | 74 | 02-8005 | 27 | 07-5001 | 53 | 21-0200 | 10 |
| 01-2110 74 03-2310 87 07-5010 53 23-0010 11 01-2115 74 03-3100 87 07-5025 53 23-0011 11 11 01-2112 74 03-3100 87 07-5025 53 23-0011 11 11 01-2122 74 03-322 85 07-6012 52 23-0200 11 10 1-2125 74 03-4221 85 07-6012 52 23-0201 11 10 1-220X 75 03-6210 88 07-6012 52 23-0000 11 10 1-221X 76 04-0008U 24 07-6048 52 23-0010S 12 01-241X 74 04-0015U 24 07-6048 52 23-0010S 12 01-241X 74 04-0015U 24 07-6096 52 23-0011S 12 01-242X 74 04-0015U 24 07-8025 50 23-0021S 12 01-242X 74 04-0050U 24 07-8025 50 23-0021S 12 01-244X 74 04-0125U 24 07-8025 50 23-0020S 12 01-270X 76 04-0500U 24 07-9075 50 23-0020S 12 01-270X 76 04-0500U 24 07-9075 50 23-0020S 12 01-2801 76 04-0100U 24 07-9175 50 23-0101S 12 01-2805 74 04-1038U 24 08-PES1322R 56 25-0051 20 01-2803 74 04-1038U 24 08-PES322R 56 25-0051 20 01-330X 77 04-1030U 24 08-PES3322R 56 25-0052 20 01-310X 77 04-1030U 24 08-PES3322R 56 25-0052 20 01-330X 77 04-1050U 24 08-PES3322R 56 25-0052 20 01-330X 77 04-1050U 24 08-PES3322R 56 25-1020 20 01-330X 77 04-1050U 24 08-PES3328R 56 30-0138 21 01-821X 81 04-2060 24 08-PES3328R 56 30-0138 21 01-821X 81 04-2060 24 08-PES3328R 56 30-0138 21 01-800X 81 04-2050 24 09-PES250F 56 30-0138 21 01-800X 81 04-2050 24 09-PES250F 56 30-0138 21 01-800X 81 04-2050 24 09-PES250F 56 60-0058 60 01-800X 86 04-3030 24 10-0850 07 51-0800 68 01-800X 86 04-3030 24 10-0850 07 51-0800 68 01-800X 86 04-3030 24 10-0850 07 51-0800 68 01-800X 86 04-3050 24 10-9850 07 51-0800 68 01-800X 86 04-3050 24 10-9850 07 51-0800 68 01-800X 86 04-3050 24 10-985 | | 74 | 02-9010 | 27 | 07-5002 | 53 | 21-1000 | 10 |
| 01-2115 | 01-2104 | 74 | 03-1300 | | | | | |
| 01-2112 | 01-2110 | 74 | | | | | | |
| 01-2122 74 03-4221 85 07-6006 52 23-0201 11 01-2155 74 03-4221 85 07-6012 52 23-1000 11 1 01-221X 76 04-0008U 24 07-6024 52 23-1001 11 1 01-221X 76 04-0008U 24 07-6086 52 23-0011S 12 01-241X 74 04-0016U 24 07-6096 52 23-0011S 12 01-241X 74 04-0060U 24 07-8025 50 23-0021S 12 01-242X 74 04-0060U 24 07-8075 50 23-0051S 12 01-243X 74 04-0060U 24 07-8075 50 23-0051S 12 01-243X 74 04-0050U 24 07-8075 50 23-0051S 12 01-243X 74 04-0050U 24 07-8075 50 23-0051S 12 01-243X 74 04-0250U 24 07-8075 50 23-0051S 12 01-243X 74 04-0250U 24 07-8075 50 23-0051S 12 01-2801 76 04-0500U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 08-PES1322R 56 25-0051 20 01-2803 74 04-1016U 24 08-PES3022R 56 25-0051 20 01-310X 77 04-1030U 24 08-PES3022R 56 25-0025 20 01-320X 77 04-1050U 24 08-PES3022R 56 25-1025 20 01-330X 77 04-1125U 24 08-PES3322S 56 25-1025 20 01-330X 77 04-1125U 24 08-PES3322S 56 25-1100 20 01-330X 77 04-1125U 24 08-PES3322S 56 25-1100 20 01-330X 77 04-1125U 24 08-PES3322M 56 30-0135 21 01-420X 89 04-1500U 24 08-PES3324M 56 30-0135 21 01-420X 89 04-1500U 24 08-PES3325M 56 30-0135 21 01-420X 89 04-1205 24 09-02250 56 41-5005 20 01-800X 81 04-2030 24 09-PES2250F 56 41-5005 60 01-80 00 01-80 | 01-2115 | 74 | | | | | | |
| 01-2155 74 03-4222 85 07-6012 52 23-1000 11 01-220X 75 03-6210 88 07-6024 52 23-1001 11 01-221X 76 04-0008U 24 07-6048 52 23-0010S 12 01-221X 74 04-0015U 24 07-6068 52 23-0011S 12 01-241X 74 04-0030U 24 07-8025 50 23-0021S 12 01-242X 74 04-0030U 24 07-8075 50 23-0021S 12 01-243X 74 04-0050U 24 07-8075 50 23-0011S 12 01-244X 74 04-0125U 24 07-8075 50 23-0011S 12 01-244X 74 04-0125U 24 07-8075 50 23-0101S 12 01-244X 74 04-0125U 24 07-8075 50 23-0101S 12 01-244X 74 04-0125U 24 07-9025 50 23-0200S 12 01-270X 76 04-0500U 24 07-9025 50 23-0200S 12 01-2801 76 04-0100U 24 07-9075 50 23-0201S 12 01-2805 74 04-1008U 24 08-PES1322R 56 25-0051 20 01-2803 74 04-1008U 24 08-PES1322R 56 25-0052 20 01-310X 77 04-1030U 24 08-PES322R 56 25-0052 20 01-310X 77 04-1030U 24 08-PES3322S 56 25-0052 20 01-330X 77 04-1050U 24 08-PES332S 56 25-1025 20 01-330X 77 04-1250U 24 08-PES332S 56 25-1020 20 01-330X 77 04-1250U 24 08-PES334SS 56 25-1202 20 01-330X 77 04-1250U 24 08-PES334SS 56 25-1202 20 01-330X 77 04-1250U 24 08-PES334SS 56 25-1202 20 01-330X 80 04-2008 24 08-PES332SD 56 30-0135A1 21 01-700X 80 04-2008 24 08-PES2550F 56 30-0135A1 21 01-700X 80 04-2008 24 09-PES2550F 56 30-0238A1 21 01-700X 80 04-2008 24 09-PES2550F 56 41-5025 26 01-833X 82 04-2105 24 09-PES2550F 56 41-5025 26 01-833X 82 04-2105 24 09-PES2550F 56 41-5025 26 01-830X 82 04-2105 24 09-PES2550F 56 41-5025 | 01-2112 | 74 | | | | | | |
| 01-220X 75 03-6210 88 07-6024 52 23-0010 11 01-221X 76 04-008U 24 07-6048 52 23-0010S 12 01-241X 74 04-0015U 24 07-6096 52 23-0011S 12 01-242X 74 04-0060U 24 07-8025 50 23-0021S 12 01-242X 74 04-0060U 24 07-8075 50 23-0051S 12 01-243X 74 04-0125U 24 07-8075 50 23-0051S 12 01-243X 74 04-0125U 24 07-8075 50 23-0051S 12 01-240X 74 04-0125U 24 07-8075 50 23-0051S 12 01-240X 74 04-0125U 24 07-8075 50 23-0101S 12 01-2801 76 04-0105U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 07-9075 50 23-0201S 12 01-2805 74 04-108U 24 08-PES1322R 56 25-0051 20 11-2803 74 04-103U 24 08-PES1322R 56 25-0051 20 11-320X 77 04-1060U 24 08-PES3022R 56 25-0052 20 11-320X 77 04-1060U 24 08-PES3022R 56 25-1025 20 11-330X 77 04-1125U 24 08-PES3022R 56 25-1020 20 11-330X 77 04-1125U 24 08-PES3322R 56 25-1100 20 11-330X 77 04-1125U 24 08-PES3322R 56 25-1100 20 11-330X 77 04-1125U 24 08-PES3322R 56 25-1100 20 11-330X 77 04-125U 24 08-PES3322R 56 30-013S 21 101-621XA 79 04-1100U 24 08-CA3022R 56 30-013S 21 101-621XA 79 04-1100U 24 08-CA3022R 56 30-013S 21 101-621XA 79 04-1100U 24 08-CA3025R 56 30-013S 21 101-800X 81 04-2030 24 08-CA3025R 56 30-013S 21 101-800X 81 04-2030 24 09-PES2250F 56 30-013S 21 101-800X 81 04-2030 24 09-PES2250F 56 30-013S 21 101-800X 81 04-2030 24 09-PES2250F 56 30-013S 21 101-800X 80 04-215 24 09-PES2250F 56 41-5000 68 101-800X 81 04-2030 24 10-0152 07 51-0500 68 101-800X 81 04-2030 24 10-0152 07 51-0500 68 101-800X 81 04-2030 24 10-0152 07 51-0500 68 101-800X 86 04-3030 24 10-0600 27 51-0600 68 01-800X 86 04-3030 24 10-0600 27 51-0600 68 01-800X 86 04-3030 24 10-0600 27 51-0600 68 01-900X 86 04-3030 24 10-0600 27 51-0600 68 01-900X 86 04-3030 24 10-9600 07 51-0600 68 01-900X 86 04-3030 24 10-9600 07 51-0600 69 00002 59 01-900X 86 04-3030 24 10-9600 07 51-0600 69 00002 59 01-900X 86 0 | 01-2122 | 74 | | | | | | |
| 01-221X 76 04-0008U 24 07-6048 52 23-0010S 12 01-241X 74 04-0015U 24 07-6096 52 23-0011S 12 01-241X 74 04-0030U 24 07-8025 50 23-0021S 12 01-243X 74 04-0050U 24 07-8075 50 23-0021S 12 01-243X 74 04-025U 24 07-8075 50 23-0011S 12 01-243X 74 04-025U 24 07-8075 50 23-0101S 12 01-240X 74 04-025U 24 07-8075 50 23-0101S 12 01-240X 74 04-025U 24 07-8075 50 23-0101S 12 01-2801 76 04-050U 24 07-9075 50 23-0200S 12 01-2801 76 04-0100U 24 07-9175 50 23-0201S 12 01-2801 76 04-0100U 24 07-9175 50 23-0201S 12 01-2801 77 04-1010U 24 07-9175 50 23-0201S 12 01-2803 74 04-1016U 24 08-PES1322R 56 25-0051 20 01-320X 77 04-1030U 24 08-PES322R 56 25-0052 20 01-310X 77 04-1030U 24 08-PES3322S 56 25-0105 20 01-320X 77 04-1060U 24 08-PES3322S 56 25-1025 20 01-330X 77 04-1125U 24 08-PES3322S 56 25-1100 20 01-330X 77 04-1125U 24 08-PES3322M 56 30-0135 21 01-202 01-330X 77 04-1125U 24 08-PES3322M 56 30-0135 21 01-202 01-330X 77 04-1250U 24 08-PES332EM 56 30-0135 21 01-202 01-330X 77 04-1250U 24 08-PES332EM 56 30-0135 21 01-700X 80 04-2008 24 08-PES3345M 56 30-0135A1 21 01-700X 80 04-2008 24 08-PES332FM 56 30-0135A1 21 01-621XA 79 04-1100U 24 08-CA3022R 56 30-0138A1 21 01-700X 80 04-2008 24 09-PES2250F 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250F 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250F 56 41-5025 26 01-830X 82 04-2105 24 09-PES2250F 56 41-5025 26 01-830X 82 04-2105 24 09-CA2250F 56 41-5005 66 01-800X 81 04-2030 24 09-CA2250F 56 41-5005 68 01-800X 82 04-2105 24 09-CA2250F 56 41-5005 60 01-800X 86 04-30330 24 10-9050 07 51-0600 68 01-800X 86 04-30330 24 10-9050 07 51-0600 68 01-900X 86 04-3015 24 10-9050 07 51-0600 68 01-900X 86 04-3015 24 10-9050 07 51-0600 68 01-900X 86 04-3015 24 10-9050 07 51-0600 68 01-900X 86 04-30330 24 10-9050 07 51-0600 68 0 | 01-2155 | 74 | | | | | | |
| 01-221X 74 04-0015U 24 07-6096 52 23-0011S 12 01-241X 74 04-0080U 24 07-8025 50 23-0021S 12 01-243X 74 04-0080U 24 07-8075 50 23-0011S 12 01-243X 74 04-0250U 24 07-8175 50 23-0011S 12 01-241X 74 04-0250U 24 07-8175 50 23-0101S 12 01-270X 76 04-050U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 07-9175 50 23-0201S 12 01-2805 74 04-1008U 24 08-PES322R 56 25-0051 20 01-2830 74 04-1015U 24 08-PES322R 56 25-0051 20 01-330X 77 04-1060U 24 08-PES3322S 56 25-1025 20 01-330X 77 04-1250U 24 08-PES3322S 56 25-1025 20 01-330X 77 04-1250U 24 08-PES332S 56 25-1020 20 01-30X 79 04-1100U 24 08-PES332S 56 25-1020 20 01-30X 80 04-2008 24 08-PES332S 56 30-0135A1 21 01-621XA 79 04-1100U 24 08-PES332S 56 30-0138A1 21 01-700X 80 04-2008 24 08-PES3250 56 30-0138A1 21 01-700X 80 04-2008 24 08-PES3250 56 30-0138A1 21 01-800X 81 04-2030 24 09-PES2250 56 30-0138A1 21 01-800X 81 04-2030 24 09-PES2250 56 30-0138A1 21 01-800X 81 04-2105 24 09-PES2250 56 30-0138A1 21 01-800X 83 04-2550 24 10-0152 07 51-0500 68 01-800X 86 04-3008 24 10-0152 07 51-0500 68 01-800X 86 04-3008 24 10-0152 07 51-0500 68 01-900X 86 04-3000 24 10-9050 07 51-0500 69 01-900X 86 | 01-220X | 75 | | | | | | |
| 01-242X 74 04-0030U 24 07-8025 50 23-0021S 12 01-243X 74 04-0105U 24 07-8075 50 23-0015IS 12 01-24X 74 04-0125U 24 07-8075 50 23-0101S 12 01-24X 74 04-0250U 24 07-9025 50 23-0200S 12 01-270X 76 04-0500U 24 07-9025 50 23-0200S 12 01-2801 76 04-0100U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 07-9175 50 23-0201S 12 01-2805 74 04-0100U 24 08-PES1322R 56 25-0051 20 01-2830 74 04-1015U 24 08-PES322R 56 25-0051 20 01-310X 77 04-1050U 24 08-PES322R 56 25-0052 20 01-310X 77 04-1050U 24 08-PES3322R 56 25-1025 20 01-330X 77 04-1125U 24 08-PES3322S 56 25-1020 20 01-330X 77 04-1125U 24 08-PES3322S 56 25-1020 20 01-330X 77 04-1250U 24 08-PES3322B 56 25-1020 20 01-330X 77 04-1125U 24 08-PES3322B 56 25-1020 20 01-330X 77 04-1125U 24 08-PES3345M 56 30-0135 21 01-420X 89 04-1500U 24 08-PES3345M 56 30-0135 21 01-420X 89 04-1500U 24 08-PES3345M 56 30-0135A1 21 01-621XA 79 04-1100U 24 08-CA3022R 56 30-0138 21 01-700X 80 04-2008 24 08-CA3045R 56 30-0138A 21 01-710X 80 04-200B 24 08-PES2250F 56 41-500 26 01-830X 81 04-2030 24 09-PES2250F 56 41-500 26 01-830X 82 04-2100 24 09-CA2250F 56 41-500 26 01-830X 82 04-2100 24 09-CA2250F 56 41-5100 26 01-830X 83 04-2500 24 10-0151 07 51-0500 68 01-830X 86 04-3015 24 10-0511 07 51-0500 68 01-830X 86 04-3015 24 10-0511 07 51-0500 68 01-830X 86 04-3015 24 10-0511 07 51-0500 68 01-830X 86 04-3015 24 10-0510 07 51-0500 68 01-930X 86 04-3050 24 10-9510 07 51-0500 69 01-930X 86 04-3050 24 10-9510 07 51-0500 69 01-930X 86 04-3500 24 10-9510 07 51-0500 69 01-930X 86 04-3500 24 10-9510 07 51-17721 69 01-950X 86 04-3500 24 10-9510 07 51-17721 69 01-950X 86 04-3500 24 10-9510 07 51-0500 69 01-950X | 01-221X | 76 | | | | | | |
| 01-243X 74 04-0060U 24 07-8075 50 23-0051S 12 01-243X 74 04-0250U 24 07-8075 50 23-0101S 12 01-244X 74 04-0250U 24 07-9025 50 23-0201S 12 01-2801 76 04-0500U 24 07-9075 50 23-0201S 12 01-2805 74 04-1008U 24 08-PES1322R 56 25-0051 20 01-2830 74 04-1015U 24 08-PES252R 56 25-0052 20 01-320X 77 04-1030U 24 08-PES3322R 56 25-0052 20 01-320X 77 04-1050U 24 08-PES3322R 56 25-1025 20 01-330X 77 04-1050U 24 08-PES3322R 56 25-1025 20 01-330X 77 04-125U 24 08-PES3322R 56 25-1025 20 01-330X 77 04-125U 24 08-PES3322R 56 25-1025 20 01-330X 77 04-125U 24 08-PES3322R 56 25-1025 20 01-320X 77 04-125U 24 08-PES3322R 56 25-1025 20 01-320X 77 04-125U 24 08-PES3322R 56 25-1025 20 01-30XA 77 04-125U 24 08-PES3322R 56 25-1025 20 01-30XA 77 04-125U 24 08-PES332R 56 25-1202 20 01-30XA 77 04-125U 24 08-PES332R 56 30-0135 21 01-420X 89 04-1500U 24 08-PES332R 56 30-0135 21 01-420X 89 04-1500U 24 08-PES332F 56 30-0135 21 01-80XA 79 04-1100U 24 08-CA3028R 56 30-0135A1 21 01-710X 80 04-2015 24 09-PES2250 56 30-038A1 21 01-710X 80 04-2015 24 09-PES2250 56 30-038A1 21 01-800X 81 04-2030 24 09-PES2250 56 41-5025 26 01-832X 82 04-215 24 09-CA2250 56 41-5025 26 01-833X 82 04-215 24 10-0151 07 51-050 68 01-870X 83 04-2250 24 10-0551 07 51-0500 68 01-870X 83 04-2250 24 10-0551 07 51-0500 68 01-870X 83 04-2250 24 10-0551 07 51-0500 68 01-870X 83 04-2250 24 10-0501 07 51-0501 68 01-900X 86 04-3030 24 10-9502 07 51-0500 68 01-970X 86 04-3030 24 10-9502 07 51-0500 68 01-970X 86 04-3030 24 10-9501 07 51-0500 68 01-970X 86 04-3030 24 10-9501 07 51-0500 68 01-970X 86 04-3030 24 10-9501 07 51-0500 68 01-970X 86 04-3050 24 10-9502 07 51-0500 68 01-970X 86 04-3050 24 10-9502 07 51-0500 68 01-970X 86 04-3030 24 10-9502 07 51-0500 68 01-970X 86 04-3030 24 10-9502 07 51-0500 68 01-970X 86 04-3050 24 10-9502 07 51-0500 68 01-970X 86 04-3050 24 10-9502 07 51-0500 68 01-970X 86 04-3050 24 10-9502 07 51-0500 68 01-970X 86 05-2721 76 12-1275 25 60-0088 60 01-970X 86 05-2721 76 12-1275 25 60-0088 60 01-970X 86 05-2721 76 12-1275 25 60-0088 60 02-0008 27 05-2793 76 12 | 01-241X | 74 | | | | | | |
| 01-244X 74 04-0125U 24 07-8175 50 23-0101S 12 01-240X 74 04-0250U 24 07-9025 50 23-0200S 12 01-2801 76 04-0500U 24 07-9075 50 23-0201S 12 01-2805 74 04-1010U 24 07-9175 50 23-1001S 12 01-2805 74 04-1010U 24 08-PES322R 56 25-0051 20 01-310X 77 04-1030U 24 08-PES322R 56 25-0052 20 01-310X 77 04-1060U 24 08-PES322R 56 25-0052 20 01-310X 77 04-1060U 24 08-PES322S 56 25-1025 20 01-320X 77 04-1060U 24 08-PES322S 56 25-1025 20 01-320X 77 04-1050U 24 08-PES322S 56 25-1025 20 01-330X 77 04-1250U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES334S 56 25-1100 20 01-320X 77 04-1250U 24 08-PES334S 56 25-1102 20 01-330X 77 04-1250U 24 08-PES334S 56 30-0135 21 01-621XA 79 04-1100U 24 08-PES334S 56 30-0135 21 01-700X 80 04-2008 24 08-PES332M 56 30-0138 21 01-710X 80 04-2015 24 09-PES2250F 56 30-0138 21 01-710X 80 04-2015 24 09-PES2250F 56 30-0138 21 01-700X 80 04-2015 24 09-PES2250F 56 41-5025 26 01-832X 82 04-2100 24 09-PES2250F 56 41-5025 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 82 04-2100 24 09-CA2250F 56 41-5025 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 83 04-2500 24 10-0501 07 51-0501 68 01-900X 86 04-3030 24 10-0502 07 51-0501 68 01-900X 86 04-3030 24 10-0502 07 51-0500 68 01-900X 86 04-3030 24 10-0502 07 51-0501 68 01-900X 86 04-3030 24 10-0502 07 51-0501 68 01-900X 86 04-3030 24 10-0850 07 51-0501 68 01-900X 86 04-3030 24 10-0850 07 51-0501 68 01-900X 86 04-3030 24 10-0850 07 51-0501 68 01-900X 86 04-3030 24 10-9850 07 51-0501 68 01-900X 86 04-3030 24 10-9850 07 51-17720 69 01-980X 86 04-3100 24 10-9850 07 51-17720 69 01-980X 86 04-3250 24 10-9850 07 51-17720 69 01-980X 86 05-2791 76 12-1299 25 60-0082 59 01-970X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2793 76 12-1399 25 60-0082 59 02-0008 27 05-2793 76 12-1399 25 60-0082 59 02-0008 27 05-2793 76 12-1399 25 60-0082 59 02-0008 27 05-2793 76 12-139 | 01-242X | 74 | | | | | | |
| 01-270X 76 04-050U 24 07-9075 50 23-020S 12 01-2801 76 04-050U 24 07-9075 50 23-0201S 12 01-2801 76 04-0100U 24 07-9075 50 23-1001S 12 01-2805 74 04-1008U 24 08-PES1322R 56 25-0051 20 01-380X 77 04-1030U 24 08-PES322R 56 25-0052 20 01-320X 77 04-1060U 24 08-PES322S 56 25-1025 20 01-320X 77 04-1060U 24 08-PES322S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES334SS 56 25-1202 20 01-320X 79 04-1100U 24 08-PES334SS 56 30-0135 21 01-420X 89 04-1500U 24 08-PES334SM 56 30-0135 21 01-621XA 79 04-1100U 24 08-PES334SM 56 30-0135A1 21 01-710X 80 04-2008 24 08-PES334SM 56 30-0138A1 21 01-710X 80 04-2015 24 09-PES2250 56 30-0138A1 21 01-800X 81 04-2030 24 09-PES2250 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250 56 41-5025 26 01-821X 81 04-2030 24 09-PES2250 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250 56 41-5100 26 01-830X 82 04-2100 24 09-CA2250 56 41-5100 26 01-830X 82 04-215 24 10-0151 07 51-0500 68 01-870X 83 04-2500 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0152 07 51-0500 68 01-870X 83 04-3008 24 10-0850 07 51-0501 68 01-990X 86 04-3015 24 10-0850 07 51-0501 68 01-990X 86 04-3008 24 10-0850 07 51-0805 70 01-990X 86 04-3008 24 10-0850 07 51-0805 70 01-990X 86 04-3060 24 10-9815 07 51-0805 70 01-990X 86 04-3060 24 10-9815 07 51-0805 70 01-990X 86 04-3250 24 10-9850 07 51-17721 69 01-990X 86 05-2721 76 12-1275 25 60-0082 59 01-990X 86 05-2721 76 12-12399 25 60-0082 59 01-990X 86 05-2721 76 12-1399 25 60-0082 59 01-990X 86 05-2721 76 12-1399 25 60-0082 59 01-990X 86 05-2721 76 12-139 | 01-243X | 74 | | | | | | |
| 01-270X 76 04-0500U 24 07-9075 50 23-0201S 12 01-2805 76 04-0100U 24 07-9175 50 23-1001S 12 01-2805 74 04-1008U 24 08-PES1322R 56 25-0051 20 01-2830 74 04-1015U 24 08-PES2522R 56 25-0052 20 01-310X 77 04-1060U 24 08-PES322R 56 25-1025 20 01-320X 77 04-1060U 24 08-PES3322S 56 25-1100 20 01-330X 77 04-1125U 24 08-PES332S 56 25-1100 20 01-330X 77 04-125U 24 08-PES3345S 56 25-1100 20 01-330X 77 04-125U 24 08-PES3345S 56 25-1202 20 01-330X 77 04-125U 24 08-PES3322M 56 30-0135 21 01-621XA 79 04-1100U 24 08-PES3322M 56 30-0138 21 01-700X 80 04-2008 24 08-CA3045R 56 30-0138 21 01-700X 80 04-2008 24 08-CA3045R 56 30-0138 21 01-710X 80 04-2015 24 09-PES250 56 30-0138 21 01-821X 81 04-2030 24 09-PES250 56 30-038A1 21 01-821X 81 04-2060 24 09-CA2250F 56 41-5025 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 82 04-2250 24 10-0152 07 51-0501 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-970X 86 04-3003 24 10-0502 07 51-0501 68 01-970X 86 04-3005 24 10-0501 07 51-0501 68 01-970X 86 04-3008 24 10-0502 07 51-0501 68 01-970X 86 04-3008 24 10-0502 07 51-0501 68 01-970X 86 04-3008 24 10-0502 07 51-0501 68 01-970X 86 04-3008 24 10-9501 07 51-0501 68 01-970X 86 04-3008 24 10-9501 07 51-0501 68 01-970X 86 04-3008 24 10-9501 07 51-0501 68 01-970X 86 04-3250 24 10-9501 07 51-0505 70 01-980X 86 04-3250 24 10-9501 07 51-0505 70 01-980X 86 04-3250 24 10-9501 07 51-17720 69 01-980X 86 04-3250 24 10-9501 07 51-0505 70 01-980X 86 04-3250 24 10-9501 07 51-0505 70 01-980X 86 04-3250 24 10-9501 07 51-0505 70 01-980X 86 04-3250 24 10-9501 07 51-17720 69 01-980X 86 04-3250 24 10-9501 07 51-17720 69 01-980X 86 04-3250 24 10-9501 07 51-17720 69 01-980X 86 05-2721 76 12-1275 25 60-0088 60 02-2791 76 12-1275 25 60-0088 60 02-2791 76 12-1275 25 60-0088 60 02-2791 76 12-1275 25 60-0088 60 02-2791 76 12-1299 25 60-0088 60 02-2791 76 12-1299 25 60-0088 60 02-2006 27 05- | 01-244X | 74 | | | | | | |
| 01-280fs 74 04-000U 24 07-9175 50 23-1001S 12 01-2830 74 04-1008U 24 08-PES1322R 56 25-0051 20 01-310X 77 04-1030U 24 08-PES322R 56 25-0052 20 01-310X 77 04-1030U 24 08-PES322R 56 25-0052 20 01-320X 77 04-1030U 24 08-PES3322R 56 25-1025 20 01-330X 77 04-1125U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES334S 56 25-1202 20 01-330XA 77 04-1250U 24 08-PES334SS 56 25-1202 20 01-330XA 77 04-1250U 24 08-PES334SS 56 25-1202 20 01-330XA 77 04-1250U 24 08-PES332SM 56 30-0135 21 01-420X 89 04-1500U 24 08-PES332SM 56 30-0135 21 01-700X 80 04-2008 24 08-CA3022R 56 30-0138 21 01-700X 80 04-2008 24 08-CA3045R 56 30-0138A1 21 01-710X 80 04-2008 24 09-PES2250 56 30-0138A1 21 01-710X 80 04-2008 24 09-PES2250 56 30-0138A1 21 01-800X 81 04-2030 24 09-PES2250 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250 56 41-5100 26 01-821X 81 04-2260 24 09-CA2250 56 41-6025 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-800X 83 04-2500 24 10-0152 07 51-0500 68 01-800X 86 04-3015 24 10-0501 07 51-0501 68 01-900X 86 04-3015 24 10-0501 07 51-0501 68 01-900X 86 04-3016 24 10-9502 07 51-0505 70 01-900X 86 04-3015 24 10-0815 07 51-0801 68 01-900X 86 04-3015 24 10-0850 07 51-0801 68 01-900X 86 04-3015 24 10-0850 07 51-0801 68 01-900X 86 04-3015 24 10-9551 07 51-0801 68 01-900X 86 04-3015 24 10-9550 07 51-0801 68 01-900X 86 04-3008 24 10-9550 07 51-0801 68 01-900X 86 04-3500 24 10-9550 07 51-0801 68 01-900X 86 04-3500 24 10-9550 07 51-0801 68 01-900X 86 05-2721 76 12-1299 25 60-0008 27 05-2793 76 12-1295 25 60-0008 60 01-900X 86 05-2791 76 12-1295 25 60-0008 60 01-900X 86 05-2791 76 12-1295 25 60-0008 60 0 | 01-270X | 76 | | | | | | |
| 01-2805 74 04-1008U 24 08-PES1322R 56 25-0051 20 01-2830 74 04-1015U 24 08-PES252R 56 25-0052 20 01-310X 77 04-1030U 24 08-PES302R 56 25-1025 20 01-320X 77 04-1060U 24 08-PES3022R 56 25-1025 20 01-320X 77 04-1125U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1125U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES332S 56 30-0135 21 01-420X 89 04-1500U 24 08-PES3345M 56 30-0135 21 01-621XA 79 04-1100U 24 08-PES3345M 56 30-0135 21 01-700X 80 04-2008 24 08-PES3345M 56 30-0138 21 01-700X 80 04-2008 24 08-CA3045R 56 30-0138 21 01-710X 80 04-2030 24 09-PES2250 56 30-0138 21 01-800X 81 04-2030 24 09-PES2250 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250F 56 41-5025 26 01-831X 82 04-2100 24 09-CA2250F 56 41-5100 26 01-830X 82 04-2105 24 09-CA2250F 56 41-5100 26 01-830X 82 04-2105 24 09-CA2250F 56 41-5100 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 83 04-2250 24 10-0151 07 51-0500 68 01-870X 83 04-2500 24 10-0502 07 51-0500 68 01-870X 83 04-2500 24 10-0502 07 51-0500 68 01-900X 86 04-3030 24 10-0502 07 51-0501 68 01-900X 86 04-3030 24 10-9502 07 51-0500 68 01-900X 86 04-3050 24 10-9502 07 51-0500 50 06 01-900X 86 04-3050 24 10-9502 07 51-0500 68 01-900X 86 04-3050 24 10-9502 07 51-0500 50 00-000 50 01-900X 86 04-3050 24 10-9502 07 51-0500 50 00-000 50 01-900X 86 05-2720 76 12-1299 25 60-0008 50 01-900X 86 05-2720 76 12-1399 25 60-0008 50 01-900X 86 05-2721 76 12-1275 25 60-0008 60 01-900X 86 05-2721 76 12-1275 25 60-0008 60 01-900X 86 05-2721 | 01-2801 | 76 | | | | | | |
| 01-2830 74 04-1015U 24 08-PES2522R 56 25-0052 20 01-310X 77 04-1030U 24 08-PES3022R 56 25-1025 20 01-320X 77 04-1050U 24 08-PES3022R 56 25-1100 20 01-330X 77 04-1125U 24 08-PES3322S 56 25-1100 20 01-330X 77 04-1250U 24 08-PES3322S 56 25-11202 20 01-330X 77 04-1250U 24 08-PES3345S 56 25-1202 20 01-330X 77 04-1250U 24 08-PES3345M 56 30-01355 21 01-420X 89 04-1500U 24 08-PES3345M 56 30-0135A1 21 01-621XA 79 04-1100U 24 08-PES3345M 56 30-0135A1 21 01-700X 80 04-2008 24 08-CA3022R 56 30-0138A1 21 01-710X 80 04-2015 24 09-PES2250 56 30-0238A1 21 01-710X 80 04-2030 24 09-PES2250F 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250F 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250F 56 41-6025 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-900X 86 04-3015 24 10-0850 07 51-0501 68 01-900X 86 04-3015 24 10-0815 07 51-0801 68 01-900X 86 04-3030 24 10-9151 07 51-0801 68 01-900X 86 04-3030 24 10-9151 07 51-0801 68 01-900X 86 04-3030 24 10-9151 07 51-0801 68 01-900X 86 04-3030 24 10-9850 07 51-0801 68 01-900X 86 04-3030 24 10-9850 07 51-0801 68 01-900X 86 04-3030 24 10-9850 07 51-0801 68 01-900X 86 04-3050 24 10-9850 07 51-17720 69 01-900X 86 04-3250 24 10-9850 07 51-17720 69 01-900X 86 05-2721 76 12-1275 25 60-0008 50 01-900X 86 05-2721 76 12-1275 25 60-0008 60 01-900X 86 05-2791 76 12-1275 25 60-0008 50 00-0025 59 02-1016 27 07-2101 54 | 01-2805 | 74 | | | | | | |
| 01-310X 77 | 01-2830 | 74 | | | | | | |
| 01-320X 77 04-1060U 24 08-PES332S 56 25-1100 20 01-330X 77 04-1125U 24 08-PES332SS 56 25-1202 20 01-330XA 77 04-1250U 24 08-PES334SS 56 25-1202 20 01-320XA 77 04-1250U 24 08-PES334SM 56 30-0135 21 01-420X 89 04-1500U 24 08-PES334SM 56 30-0135A1 21 01-621XA 79 04-1100U 24 08-CA3022R 56 30-0138 21 01-700X 80 04-2008 24 08-CA3022R 56 30-0138 21 01-710X 80 04-2015 24 09-PES2250 56 30-0238A1 21 01-710X 80 04-2015 24 09-PES2250F 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250F 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250F 56 41-5025 26 01-821X 81 04-2100 24 09-CA2250F 56 41-6025 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0151 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0500 68 01-970X 86 04-3030 24 10-0501 07 51-0501 68 01-910X 86 04-3030 24 10-0815 07 51-0800 68 01-910X 84 04-3008 24 10-0815 07 51-0800 68 01-910X 84 04-3030 24 10-9850 07 51-0801 68 01-910X 84 04-3030 24 10-9151 07 51-0801 68 01-910X 84 04-3030 24 10-9501 07 51-0801 68 01-910X 84 04-3030 24 10-9501 07 51-0801 68 01-910X 84 04-3030 24 10-9850 07 51-0801 68 01-910X 84 04-3030 24 10-9501 07 51-0801 68 01-910X 84 04-3030 24 10-9501 07 51-0805 70 01-920X 86 04-3125 24 10-9151 07 51-0805 60 01-980X 86 04-3125 24 10-9501 07 51-17720 69 01-930X 86 04-3125 24 10-9501 07 51-0805 70 01-920X 86 04-3125 24 10-9501 07 51-17721 69 01-980X 86 04-3250 24 10-9502 07 51-17720 69 01-980X 86 05-2721 76 12-1275 25 60-0008 60 01-980X 86 05-2721 76 12-1275 25 60-0082 59 01-970X 86 05-2721 76 12-1275 25 60-0082 60 01-980X 86 05-2721 76 12-1275 25 60-0082 60 01-980X 86 05-2721 76 12-1279 25 60-0082 69 02-0088 27 05-2793 76 12-1399 25 60-0082 59 02-0088 27 05-2793 76 12-1399 25 60-0082 59 02-0088 27 05-2793 76 12-1399 25 60-0082 59 02-1016 27 05-2793 76 12-1399 25 60-0082 59 02-1016 27 05-2793 76 12-1399 25 60-0082 59 02-1016 27 05-2793 76 12-1399 25 60-0082 59 02-1016 27 05-2793 76 12-1399 25 60-0108 60 02-1016 27 05-2793 76 12-1399 25 60-0108 60 02-1016 27 05-2793 76 12-1399 25 60-0108 60 02-1016 27 05-2793 76 12-1299 25 60-0108 60 02-101 | 01-310X | 77 | | | | | | |
| 01-330X | 01-320X | 77 | | | | | | |
| 01-330XA 77 | 01-330X | 77 | | | | | | |
| 01-420X 89 04-1500U 24 08-PES3345M 56 30-0135A1 21 01-621XA 79 04-1100U 24 08-CA3022R 56 30-0138 21 01-700X 80 04-2008 24 08-CA3045R 56 30-0138A1 21 01-710X 80 04-2015 24 09-PES2250 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250 56 30-0238A1 21 01-800X 81 04-2060 24 09-PES2250F 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250 56 41-5100 26 01-830X 82 04-2100 24 09-CA2250F 56 41-6100 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 83 04-2500 24 10-0151 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-870X 83 04-3008 24 10-0502 07 51-0505 70 01-900X 86 04-3015 24 10-0815 07 51-0800 68 01-910X 84 04-3030 24 10-0850 07 51-0800 68 01-910X 84 04-3030 24 10-0850 07 51-0800 68 01-910X 84 04-3030 24 10-9151 07 51-0801 68 01-930X 85 04-3125 24 10-9151 07 51-0801 68 01-930X 85 04-3125 24 10-9501 07 51-0805 70 01-930X 86 04-3015 24 10-9501 07 51-0800 68 01-930X 85 04-3125 24 10-9501 07 51-0805 70 01-930X 86 04-3315 24 10-9501 07 51-0805 68 01-930X 85 04-3125 24 10-9501 07 51-0805 70 01-930X 85 04-3125 24 10-9501 07 51-17720 69 01-930X 86 04-3250 24 10-9501 07 51-17720 69 01-930X 86 04-3250 24 10-9501 07 51-17725 70 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 12-1279 25 60-0082 59 01-970X 86 05-2721 76 12-1279 25 60-0082 59 01-970X 86 05-2721 76 12-1279 25 60-0082 59 01-970X 86 05-2721 76 12-1279 25 60-0082 59 01-990X 86 05-2721 76 12-1299 25 60-0082 59 02-0008 27 05-2793 76 12-1299 25 60-0082 59 02-0008 27 05-2793 76 12-1299 25 60-0082 59 02-0008 27 05-2793 76 12-1299 25 60-0082 59 02-0008 27 05-2793 76 12-1299 25 60-0102 59 02-1016 27 05-2793 76 12-1275 25 60-0102 59 02-1016 27 05-2793 76 12-1275 25 60-0108 60 02-1022 27 05-2793 76 12-1270 55 60-0188 60 02-1016 27 07-2101 54 15-1070 55 60-0188 60 02-1016 27 07-2101 54 15-1070 55 60-0188 60 02-1016 27 07-2102 54 15-1100 55 60-0188 60 02-1016 27 07-2104 54 20-0010 10 60-0208 50 02-1016 27 07-2108 54 20-0010 10 60-0208 59 02-1016 27 07-2108 54 20-0010 10 60-0208 59 02-1016 27 07-2108 54 20-0010 10 60-0208 50 02-1016 27 07-2108 54 20-00 | 01-330XA | 77 | | | | | | |
| 01-621XA 79 04-1100U 24 08-CA3022R 56 30-0138 21 01-700X 80 04-2008 24 08-CA3045R 56 30-0138A1 21 01-710X 80 04-2015 24 09-PES2250 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250F 56 41-5005 26 01-821X 81 04-2060 24 09-CA2250F 56 41-5100 26 01-830X 82 04-2100 24 09-CA2250F 56 41-5100 26 01-830X 82 04-2125 24 10-0151 07 41-6100 26 01-830X 82 04-2125 24 10-0152 07 51-0500 68 01-860X 83 04-2500 24 10-0501 07 51-0501 68 01-90X 86 04-3015 24 10-0502 07 51-0505 70 01-90X <td>01-420X</td> <td>89</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | 01-420X | 89 | | | | | | |
| 01-700X 80 04-2008 24 08-CA3045R 56 30-0138A1 21 01-710X 80 04-2015 24 09-PES2250 56 30-0238A1 21 01-800X 81 04-2060 24 09-PES2250F 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250F 56 41-5100 26 01-830X 82 04-2100 24 09-CA2250F 56 41-6100 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0501 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-870X 83 04-3008 24 10-0502 07 51-0505 70 01-90X 86 04-3015 24 10-0815 07 51-0505 70 01-90X | 01-621XA | 79 | | 24 | | | | |
| 01-710X 80 04-2015 24 09-PES2250 56 30-0238A1 21 01-800X 81 04-2030 24 09-PES2250F 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250F 56 41-5100 26 01-830X 82 04-2105 24 09-CA2250F 56 41-6100 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-900X 86 04-3008 24 10-0502 07 51-0501 68 01-90X 86 04-3030 24 10-0850 07 51-0801 68 01-940X 86 04-3103 24 10-9151 07 51-0806 70 01-920X | | 80 | | 24 | | | | |
| 01-800X 81 04-2030 24 09-PES2250F 56 41-5025 26 01-821X 81 04-2060 24 09-CA2250F 56 41-5100 26 01-830X 82 04-2100 24 09-CA2250F 56 41-6025 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-900X 86 04-3015 24 10-0815 07 51-0805 70 01-90X 86 04-3015 24 10-0815 07 51-0806 68 01-940X 86 04-3015 24 10-0815 07 51-0805 70 01-920X 84 04-3060 24 10-9151 07 51-0805 70 01-930X <td< td=""><td></td><td>80</td><td>04-2015</td><td>24</td><td>09-PES2250</td><td>56</td><td>30-0238A1</td><td>21</td></td<> | | 80 | 04-2015 | 24 | 09-PES2250 | 56 | 30-0238A1 | 21 |
| 01-821X 81 04-2060 24 09-CA2250 56 41-5100 26 01-830X 82 04-2100 24 09-CA2250F 56 41-6025 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0505 70 01-900X 86 04-3008 24 10-0502 07 51-0505 70 01-90X 86 04-3015 24 10-0815 07 51-0800 68 01-9A0X 86 04-3030 24 10-0850 07 51-0801 68 01-9A0X 84 04-3100 24 10-9152 07 51-07005 70 01-920X 84 04-3125 24 10-9501 07 51-17720 69 01-930X 85 | | | 04-2030 | 24 | 09-PES2250F | 56 | 41-5025 | 26 |
| 01-830X 82 04-2100 24 09-CA2250F 56 41-6025 26 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-900X 86 04-3008 24 10-0502 07 51-0505 70 01-900X 86 04-3015 24 10-0815 07 51-0800 68 01-9A0X 86 04-3030 24 10-0850 07 51-0800 68 01-9A0X 84 04-3060 24 10-9151 07 51-0805 70 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86< | | | 04-2060 | 24 | 09-CA2250 | 56 | 41-5100 | 26 |
| 01-832X 82 04-2125 24 10-0151 07 41-6100 26 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-900X 86 04-3008 24 10-0502 07 51-0505 70 01-90X 86 04-3015 24 10-0815 07 51-0800 68 01-9A0X 86 04-3030 24 10-0850 07 51-0801 68 01-910X 84 04-3060 24 10-9151 07 51-0805 70 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86 04-3500 24 10-9815 07 60-0002 59 01-960X 86 | | | 04-2100 | 24 | 09-CA2250F | 56 | 41-6025 | 26 |
| 01-860X 83 04-2250 24 10-0152 07 51-0500 68 01-870X 83 04-2500 24 10-0501 07 51-0501 68 01-900X 86 04-3008 24 10-0502 07 51-0505 70 01-9A0X 86 04-3030 24 10-0815 07 51-0801 68 01-910X 84 04-3060 24 10-9151 07 51-0805 70 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-930X 86 04-3250 24 10-9501 07 51-17725 70 01-960X 86 04-3250 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0002 59 01-980X 86 <td></td> <td></td> <td>04-2125</td> <td>24</td> <td>10-0151</td> <td>07</td> <td>41-6100</td> <td>26</td> | | | 04-2125 | 24 | 10-0151 | 07 | 41-6100 | 26 |
| 01-870X 83 04-2500 24 10-0501 07 51-0505 70 01-900X 86 04-3008 24 10-0502 07 51-0505 70 01-9A0X 86 04-3015 24 10-0815 07 51-0800 68 01-9A0X 86 04-3030 24 10-0850 07 51-0801 68 01-910X 84 04-3060 24 10-9151 07 51-0805 70 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86 04-3250 24 10-9502 07 51-17721 69 01-960X 86 04-3250 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0008 60 01-980X 86 <td></td> <td></td> <td></td> <td></td> <td>10-0152</td> <td></td> <td>51-0500</td> <td></td> | | | | | 10-0152 | | 51-0500 | |
| 01-900X 86 04-3015 24 10-0815 07 51-0800 68 01-9A0X 86 04-3030 24 10-0850 07 51-0801 68 01-910X 84 04-3060 24 10-9151 07 51-0805 70 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86 04-3250 24 10-9502 07 51-17725 70 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0002 59 01-980X 86 05-2721 76 12-1275 25 60-0088 60 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | | | | | | | | |
| 01-9A0X 86 04-3030 24 10-0815 07 51-0800 68 01-910X 84 04-3060 24 10-0850 07 51-0801 68 01-920X 84 04-3060 24 10-9151 07 51-0805 70 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86 04-3250 24 10-9502 07 51-17725 70 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0008 60 01-980X 86 05-2721 76 12-1275 25 60-0056 61 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | | | | | | | | |
| 01-910X 84 04-3060 24 10-9151 07 51-0805 70 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86 04-3250 24 10-9502 07 51-17725 70 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0002 59 01-980X 86 05-2721 76 12-1275 25 60-0008 60 01-980X 86 05-2721 76 12-1275 25 60-0082 59 02-0008 27 05-2791 76 12-1299 25 60-0082 59 02-0096 27 05-2793 76 12-1300 25 60-0102 59 02-1016 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | | | | | | | | |
| 01-920X 84 04-3100 24 10-9152 07 51-17720 69 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86 04-3250 24 10-9502 07 51-17725 70 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0002 59 01-980X 86 05-2721 76 12-1275 25 60-0008 60 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2792 76 12-1300 25 60-0088 60 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0108 60 02-1022 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | | | | | | | | |
| 01-930X 85 04-3125 24 10-9501 07 51-17721 69 01-950X 86 04-3250 24 10-9502 07 51-17725 70 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0008 60 01-980X 86 05-2721 76 12-1275 25 60-0056 61 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2792 76 12-1300 25 60-0082 59 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0108 60 02-1022 27 05-2798 76 15-1040 55 60-0186 61 02-1096 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | | | | | | | | |
| 01-950X 86 04-3250 24 10-9502 07 51-17725 70 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0008 60 01-980X 86 05-2721 76 12-1275 25 60-0056 61 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2792 76 12-1300 25 60-0088 60 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0102 59 02-1022 27 05-2798 76 15-1040 55 60-0188 60 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 | | | | | | | | |
| 01-960X 86 04-3500 24 10-9815 07 60-0002 59 01-970X 86 05-2720 76 10-9850 07 60-0008 60 01-980X 86 05-2721 76 12-1275 25 60-0056 61 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2792 76 12-1300 25 60-0088 60 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0102 59 02-1022 27 05-2798 76 15-1040 55 60-0108 60 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | | | | | | | | |
| 01-970X 86 05-2720 76 10-9850 07 60-0008 60 01-980X 86 05-2721 76 12-1275 25 60-0056 61 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2792 76 12-1300 25 60-0088 60 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0102 59 02-1022 27 05-2798 76 15-1040 55 60-0108 60 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> | | | | | | | | |
| 01-980X 86 05-2721 76 12-1275 25 60-0056 61 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2792 76 12-1300 25 60-0088 60 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0108 60 02-1022 27 05-2798 76 15-1040 55 60-0156 61 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0182 59 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 01-990X 86 05-2791 76 12-1299 25 60-0082 59 02-0008 27 05-2792 76 12-1300 25 60-0088 60 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0108 60 02-1022 27 05-2798 76 15-1040 55 60-0156 61 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0182 59 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-0008 27 05-2792 76 12-1300 25 60-0088 60 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0108 60 02-1022 27 05-2798 76 15-1040 55 60-0156 61 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-0096 27 05-2793 76 12-1399 25 60-0102 59 02-1016 27 05-2797 76 12-2275 25 60-0108 60 02-1022 27 05-2798 76 15-1040 55 60-0156 61 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-1016 27 05-2797 76 12-2275 25 60-0108 60 02-1022 27 05-2798 76 15-1040 55 60-0156 61 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-1022 27 05-2798 76 15-1040 55 60-0156 61 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-1096 27 07-2101 54 15-1070 55 60-0182 59 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-1196 27 07-2102 54 15-1100 55 60-0188 60 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-2010-50 27 07-2104 54 20-0010 10 60-0202 59 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02-2196 27 07-2108 54 20-0200 10 60-0208 60 | | | | | | | | |
| 02 2100 | | | | | | | | |
| | | | | | | | | |

Catalog Number Index

| Cat # | Page# | Cat # | Page# | Cat # | Page# | Cat # | Page# |
|-----------|----------|--------------------|----------|---------|----------|--------------------------|----------|
| 60-0282 | 59 | 81-0056 | 31 | 89-910X | 35 | 97–6115 | 29 |
| 60-0288 | 60 | 81-0153 | 31 | 89-915X | 35 | 98-0113 | 40 |
| 60-0302 | 59 | 81-0154 | 31 | 89-920X | 35 | 98-0200 | 40 |
| 60-0308 | 60 | 81-0203 | 31 | 89-921X | 35 | 98-0213 | 40 |
| 60-0356 | 61 | 81-0204 | 31 | 89-922X | 35 | 98-0313 | 40 |
| 60-0382 | 59 | 81-0205 | 31 | 89-923X | 35 | 98-1213 | 40 |
| 60-0388 | 60 | 81-7053 | 32 | 89-950X | 35 | 98-4127 | 36 |
| 60-0402 | 59 | 81-7054 | 32 | 90-1200 | 38 | 98-4227 | 36 |
| 60-0408 | 60 | 81–7055 | 32 | 90-1281 | 38 | 98-5227 | 37 |
| 60-0456 | 61 | 81-7056 | 32 | 90-1536 | 09 | 99-2215 | 42 |
| 60-0482 | 59 | 81-7153 | 32 | 90-1550 | 09 | 99-2216 | 42 |
| 60-0488 | 60 | 81-7154 | 32 | 90-2200 | 38 | 99-2220 | 42 |
| 61-0010 | 28 | 81-7203 | 32 | 90-2281 | 38 | 99-2224 | 42 |
| 61-0020 | 28 | 81-7204 | 32 | 90-2300 | 38 | 99-2225 | 42 |
| 61-0040 | 28 | 81-7205 | 32 | 90-2381 | 38 | 99-2230 | 42 |
| 61-0041 | 28 | 81-8204 | 32 | 90-5016 | 09 | 99-5216 | 43 |
| 61-0051 | 28 | 88-005X | 34 | 90-5025 | 09 | 99-5220 | 43 |
| 61-0060 | 28 | 88-010X | 34 | 90-5200 | 38 | 99-5224 | 43 |
| 61-0080 | 28 | 88-015X | 34 | 90-5281 | 38 | 99-5225 | 43 |
| 61-0120 | 28 | 88-020X | 34 | 90-3281 | 09 | 99-62XX | 43 |
| 61-0152 | 28 | 88-020X 88-021X | 34 | 90-8009 | 39 | 99-7216 | 42 |
| 61-0132 | 28 | 88-050X | 34 | 90-8100 | 39 | 99-7220 | 42 |
| 61-0210 | 28 | 88-321X | 33 | 90-8281 | 39 | 99-7224 | 42 |
| | | | | | | CKS-1008 | 44 |
| 61-0304 | 28 | 88-320X | 33 33 | 90-8381 | 39 39 | CKS-1520 | 45 |
| 65-0001 | 65 65 | 88-350X | | 90-9009 | | CKS-152X | 45 |
| 65-0002 | 65 65 | 88-320XS | 33 | 90-9025 | 41 41 | CKS-0101 | 44 |
| 65-0010 | 66 | 88-350XS | 33 | 90-9081 | 41 | CKS-1301 | 44 |
| 65-1001 | 66 | 88-403X | 36 | 90-9100 | 41 | AB23-0020S | 20 |
| 65-1010 | | 88-405X | 36 | 90-9200 | | AB23-0050S | 20 |
| 66-1501 | 64 | 88-408X | 36 36 | 90-9210 | 40 | AB23-0250S | 20 |
| 66-1515 | 64 | 88-410X | 36 | 90-9250 | 40 | AT23-0020S | 20 |
| 66-1560 | 64 | 88-514X | 37 | 90-9381 | 40 | AT23-0020CS | 20 |
| 66-1701 | 64 | 88-600X | 33 | 93-0050 | 23 | AT23-0050S | 20 |
| 70-1180 | 57 | 88-605X | 34 | 93-0100 | 23 | AT23-0050CS | 20 |
| 70-1250 | 57 | 88-610X | 34 | 93-0250 | 23 | AT23-0200S | 20 |
| 70-2180 | 57 | 88-615X | 34 | 93-0400 | 23 | AT23-0200CS | 20 |
| 7101 | 26 | 88-620X | 33 | 93-0800 | 23 | AT23-1000S | 20 |
| 7105 | 26 | 88-620XS | 33 | 93-1000 | 23 | AT23-1000CS | 20 |
| 80-0015 | 80 | 88-621X | 33 | 96-2216 | 42 | AH23-0050S | 20 |
| 80-0020 | 08 | 88-650X | 33 | 96-2220 | 42 | AH23-0050CS | 20 |
| 80-0500 | 08 | 88-650XS | 33 | 96-2224 | 42 | AH23-0300S | 20 |
| 80-1500 | 08 | 88-720X | 34 | 96-62XX | 43 | AH23-0300CS | 20 |
| 80-1500-S | 08 | 88-721X | 34 | 96-5216 | 43 | AH23-1000S | 20 |
| 80-1501 | 08 | 88-750X | 34 | 96-5220 | 43 | AH23-1000CS | 20 |
| 80-1502 | 08 | 88-905X | 34 | 96-7216 | 42 | AA21-0030S | 20 |
| 80-1503 | 08 | 88-910X | 34 | 96-7220 | 42 | AA21 00303 AA23-0250S | 20 |
| 81-0000 | 31 | 88-915X | 34 | 96-7224 | 42 | R23-0020 | 19 |
| 81-0001 | 31 | 88-920X | 34 | 97–1111 | 29 | R23-0020 R23-0021 | 19 |
| 81-0002 | 31 | 88-921X | 34 | 97–1112 | 29 | | 19 |
| 81-0003 | 31 | 88-950X | 34 | 97–1113 | 29 | R23-0200 R23-0201 | 19 |
| 81-0004 | 31 | 89-423X | 36 | 97–1114 | 29 | | |
| 81-0005 | 31 | 89-425X | 36 | 97–1115 | 29 | R23-1000 R23-1001 | 19 19 |
| 81-0006 | 31 | 89-428X | 36 | 97–6111 | 29 | N23-1001 | 19 |
| 81-0053 | 31 | 89-421X | 36 | 97–6112 | 29 | | |
| 81-0054 | 31 | 89-514X | 37 | 97–6113 | 29 | | |
| 81-0055 | 31 | 89-905X | 35 | 97–6114 | 29 | | |



The General Laboratory Supplies section includes many consumables commonly used in basic laboratory and medical procedures. From centrifuge tubes and test tubes to disposable gloves and microscope slides, many of these general laboratory supplies may have applications in multiple life science fields.

Centrifuge Tubes





>> Features

- Tubes: made of strong, highly transparent medical grade polypropylene
- Flat-top caps: made of medical grade polyethylene
- Plug-seal screw on caps: made of medical grade polypropylene
- 15ml and 50ml conical tubes available;
- Large white writing area for easy marking
- Double threaded design reduces cross-threading; caps are easily opened and closed
- Flat-top caps can be written on
- Gamma radiation sterilized, non-pyrogenic, RNase&DNase free
- Available in racked package (sterile) or bulk package (sterile & non-sterile)

Specifications


- Working temperature: stable from -20°C to 121°C (except flat-top cap)
- 15ml and 50ml tubes can withstand centrifugation of 8,400 and 9,400 G respectively
- Sealing test: complies with the International Air Transport Association (IATA) Safety Standards





Black graduation marks

| Cat # | Volume | Caps | Sterile | Package | Unit | Total |
|---------|--------|-----------|---------|---------|--------------------------|-----------|
| 10-9151 | 15ml | Flat Top | Yes | Rack | 50 /Rack, 10 Racks /Case | 500 /Case |
| 10-9152 | 15ml | Flat Top | Yes | Bulk | 25 /Bag, 20 Bags /Case | 500 /Case |
| 10-9501 | 50ml | Flat Top | Yes | Rack | 25 /Rack, 20 Racks /Case | 500 /Case |
| 10-9502 | 50ml | Flat Top | Yes | Bulk | 25 /Bag, 20 Bags /Case | 500 /Case |
| 10-9815 | 15ml | Flat Top | No | Bulk | 500 /Case | 500 /Case |
| 10-9850 | 50ml | Flat Top | No | Bulk | 500 /Case | 500 /Case |
| 10-0151 | 15ml | Plug Seal | Yes | Rack | 50 /Rack, 10 Racks /Case | 500 /Case |
| 10-0152 | 15ml | Plug Seal | Yes | Bulk | 25 /Bag, 20 Bags /Case | 500 /Case |
| 10-0501 | 50ml | Plug Seal | Yes | Rack | 25 /Rack, 20 Racks /Case | 500 /Case |
| 10-0502 | 50ml | Plug Seal | Yes | Bulk | 25 /Bag, 20 Bags /Case | 500 /Case |
| 10-0815 | 15ml | Plug Seal | No | Bulk | 500 /Case | 500 /Case |
| 10-0850 | 50ml | Plug Seal | No | Bulk | 500 /Case | 500 /Case |



Dimensions = mm

Microcentrifuge Tubes





>> Features


- Made of high quality medical grade polypropylene
- Transparent tube walls permit easy viewing of contents
- Beveled lip for convenient one-handed operation
- Writing area and attached flat cap for labeling

Specifications

- Working temperature: stable from -20°C to 121°C
- Can withstand high speed centrifugation of 20,000 G



80-0015




Dimensions = mm

| Cat # | Volume | Color | Sterile | RNase & DNase free | Unit | Total |
|-----------|--------|--------|---------|--------------------|------------------------------|------------|
| 80-0500 | 0.5ml | Clear | No | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |
| 80-0015 | 1.5ml | Clear | No | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |
| 80-1500 | 1.5ml | Clear | No | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |
| 80-1500-S | 1.5ml | Clear | Yes | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |
| 80-0020 | 2.0ml | Clear | No | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |
| 80-1501* | 1.5ml | Red | No | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |
| 80-1502* | 1.5ml | Yellow | No | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |
| 80-1503* | 1.5ml | Blue | No | Yes | 500 /Bag/Pack, 10 Packs/Case | 5000 /Case |

^{*} Call for details

Centrifuge Tube Racks





Microcentrifuge Tube Racks


- Made of high strength plastic, chemically resistant to alcohols and mild organic solvents
- Writing area on side of rack for easy marking
- Handles on each end for secure handling
- 5×16 array, 80 wells, molded with a numeric grid for easy tube identification
- Fit for 1.5ml, 2.0ml microcentrifuge tubes and screw cap tubes


PP Centrifuge Tube Racks

- Made of high quality medical grade polypropylene
- 50-well racks hold 15ml centrifuge tubes
- 25-well racks hold 50ml centrifuge tubes
- Sturdy design keeps tubes upright

Cardboard Centrifuge Tube Racks

- Made of durable cardboard with superior white coating, chemically resistant to alcohols and mild organic solvents
- 36-well racks for 15ml centrifuge tubes, 16-well racks for 50ml centrifuge tubes
- Removable lid features a lined writing area





Dimensions = mm

| Cat # | Well | Color | Sterile | Unit | Total |
|---------|---------|-----------------|---------|----------------------|-----------|
| 90-8009 | 80-well | Assorted Colors | No | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-1536 | 36-well | White | No | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-5016 | 16-well | White | No | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-1550 | 50-well | Blue | No | 1 /Bag, 20 Bags/Case | 20 /Case |
| 90-5025 | 25-well | Blue | No | 1 /Bag, 20 Bags/Case | 20 /Case |




Cardboard Centrifuge Tube Rack Working temperature: from -196 °C to 121°C

Pipet Tips





- Made of polypropylene, autoclavable
- Available in 10µl, 200µl and 1000µl sizes, suitable for standard manual and electronic pipettes
- 10µl and 200µl tips are graduated, 200µl tips have a beveled point design that reduces sample retention Extended 1000µl tips facilitate pipetting from deep or narrow containers
- RNase & DNase free, non-pyrogenic
- Available in racked package (sterile) and bulk package (non-sterile)



Dimensions = mm

| Cat # | Volume | Color | Sterile | Unit | Total |
|---------|-------------------|---------|---------|--|-------------|
| 20-0010 | 10µl | Natural | No | 1,000 /Bag, 2 Bags /Pack, 5 Packs/Case | 10000 /Case |
| 21-0010 | 10µl | Natural | Yes | 96 /Rack, 10 Racks/Pack, 10 Packs/Case | 9600 /Case |
| 20-0200 | 200µl | Yellow | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 21-0200 | 200µl | Yellow | Yes | 96 /Rack, 10 Racks/Pack, 10 Packs/Case | 9600 /Case |
| 20-1000 | 1000µl | Blue | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 21-1000 | 1000µl | Blue | Yes | 100 /Rack, 6 Racks/Pack, 10 Packs/Case | 6000 /Case |
| 20-1250 | 1000µl Extended | Natural | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 21-1250 | 1000µl Extended | Natural | Yes | 96 /Rack, 6 Racks/Pack, 10 Packs/Case | 5760 /Case |
| 20-1251 | 1000µl Extra-Long | Blue | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |

Filter Tips






>> Features

- Tips are made of medical grade polypropylene, and filters are made of highly hydrophobic UHMWPE
- Filters prevent cross-contamination and damage to samples, ensure accurate and precise sample aspiration, and protect the user by preventing sample liquid or vapor from entering the pipette body
- Standard filter tips and low retention filter tips available in all sizes
- Low retention treated tips undergo a coating-free treatment and will not contaminate samples
- All filter tips are autoclavable
- Available in racked package (sterile) and bulk package (non-sterile)
- RNase & DNase free, non-pyrogenic



Small and uniform filter pore size ensures the accuracy and precision of sample aspiration



10µl Dimensions=mm

| Cat # | Volume | Color | Sterile | Low Retention | Unit | Total |
|---------|--------|---------|---------|---------------|--|-------------|
| 22-0010 | 10µl | Natural | No | No | 1,000 /Bag, 2 Bags/Pack, 5 Packs /Case | 10000 /Case |
| 22-0011 | 10µl | Natural | No | Yes | 1,000 /Bag, 2 Bags/Pack, 5 Packs /Case | 10000 /Case |
| 23-0010 | 10µl | Natural | Yes | No | 96 /Rack, 10 Racks/Pack, 10 Packs/Case | 9600 /Case |
| 23-0011 | 10µl | Natural | Yes | Yes | 96 /Rack, 10 Racks/Pack, 10 Packs/Case | 9600 /Case |
| 22-0200 | 200µl | Natural | No | No | 1,000 /Bag/Pack, 10 Packs /Case | 10000 /Case |
| 22-0201 | 200µl | Natural | No | Yes | 1,000 /Bag/Pack, 10 Packs /Case | 10000 /Case |
| 23-0200 | 200µl | Natural | Yes | No | 96 /Rack, 10 Racks/Pack, 10 Packs/Case | 9600 /Case |
| 23-0201 | 200µl | Natural | Yes | Yes | 96 /Rack, 10 Racks/Pack, 10 Packs/Case | 9600 /Case |
| 22-1000 | 1000µl | Natural | No | No | 1,000 /Bag/Pack, 10 Packs /Case | 10000 /Case |
| 22-1001 | 1000µl | Natural | No | Yes | 1,000 /Bag/Pack, 10 Packs /Case | 10000 /Case |
| 23-1000 | 1000µl | Natural | Yes | No | 96 /Rack, 6 Racks/Pack, 10 Packs/Case | 5760 /Case |
| 23-1001 | 1000µl | Natural | Yes | Yes | 96 /Rack, 6 Racks/Pack, 10 Packs/Case | 5760 /Case |


Filter Tips-S Series






>> Features


- Made in Bioclean room
- Extremely straight low retention tips ensure accurate and precise sample
- Extra-long tips: 10μl with max volume 20μl, 200μl with max volume 300μl, and 1000μl with max volume 1250μl
- Racked filter tips, sterile
- DNase & RNase free, non-pyrogenic





Dimensions = mm





| Cat # | Volume | Color | Sterile | Unit | Total |
|----------|-----------------|---------|---------|-------------------------------|------------|
| 23-00105 | 10µl Extended | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |
| 23-00115 | 10µl Extended | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |
| 23-00215 | 20µl | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |
| 23-00515 | 50µl | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |
| 23-01015 | 100µl | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |
| 23-02005 | 200µl Extended | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |
| 23-02015 | 200µl Extended | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |
| 23-10015 | 1000µl Extended | Natural | Yes | 96 Pieces/Rack, 50 Racks/Case | 4800 /Case |













| | | | | • | | | | | | | | | | | | | | | | | | | | | | | |
|----------------|----------------------|-------------|----------------|-----------|--------------|-----------|---------------------------|--|-------------|-------------|----------------|----------------|---------------|-------------------|-------------|---------------------------|---------------|------------------------|----------------|-------------------|--------------|-----------------------------|---------------|-------------|----------------|----------------------|-------------------|
| Pipet | te Brand | Ep | pend | orf R | Resea | arch | Plus | G | ilsor | n Pi | petn | nan(f | PG) |) | G | ilson | Pipe | etma | n(PI | 1) | (| Gils | on | Pip | etma | an(PL | _) |
| Mode | I / Range | 0.1-2.5ul | 0.5-10ul | 2-20ul | 10-100ul | 20-200ul | 100-1000ul | P2, 0.2-2ul | P10, 1-10ul | P20, 2-20ul | P100, 10-100ul | P200 20-20011 | 200, 20 20001 | P1000, 100-1000ul | P2, 0.2-2ul | P10, 1-10ul | P20, 2-20ul | P100, 10-100ul | P200, 20-200ul | P1000, 100-1000ul | P2, 0.2-2ul | D10 1-1011 | PIO, I-IOUI | P20, 2-20ul | P100, 10-100ul | P200, 20-200ul | P1000, 100-1000ul |
| Volume | Product No. | | | | | | | | | | Ĺ | Jnive | ersa | al Fit | Filt | er Ti | ips | | | | | | | | | | |
| 0.5–10ul | 23-0010S | $\sqrt{}$ | $\sqrt{}$ | | | | | | $\sqrt{}$ | | | | | | $\sqrt{}$ | $^{\sim}$ | | | | | $\sqrt{}$ | | | | | | |
| 1–20ul | 23-0011S 23-0020S | <i>√</i> | $\sqrt{}$ | √ | | | | $\sqrt{}$ | $\sqrt{}$ | √ | | | | | $\sqrt{}$ | $\sqrt{}$ | √ | | | | <i>√</i> | | | √ | | | |
| 1-20ui | 23-0021S | | | $\sqrt{}$ | L | | | | | | | | _ | _ | | | $\sqrt{}$ | | _ | | | ┸ | | $\sqrt{}$ | | | |
| 1–50ul | 23–0050S 23–0051S | | | | | | | | | | | - √ - √ | + | + | | | | $\sqrt{}$ | $\frac{}{}$ | | | + | \dashv | | | $\sqrt{}$ | |
| 1–100ul | 23-0100S | | | | $\sqrt{}$ | | | | | | | - √ | | | | | | $\sqrt{}$ | $\sqrt{}$ | | | | | | $\sqrt{}$ | $\sqrt{}$ | |
| 1-10001 | 23-0101S | | | | $\sqrt{}$ | | | | | | | - √ | 1 | | | | | $\sqrt{}$ | $\sqrt{}$ | | | | | | $\sqrt{}$ | $\sqrt{}$ | |
| 1-300ul | 23-0200S | | _ | | | \ | | | | | | | 1 | | | | | | $\sqrt{}$ | | | ┸ | | | | $\sqrt{}$ | |
| | 23-0201S | | | | | √ | _ | | | ┡ | _ | | | _ | | | | | $\sqrt{}$ | _ | | ╄ | | | | $\sqrt{}$ | _ |
| 100- | 23-1000S | | - | | | | √ | | | ┡ | + | + | | √ | | | | | | √ | | ╄ | | | | | √ |
| 1000ul | 23-1001S | | | | | | - √ | | | | | | | √ | | | | | | √ | | | | | | | √ |
| Pinet | te Brand | 1 | | nerm | no Fir | nnnir | ette | (F1) |) | | | Ther | rmo | o Fin | npir | ette | (F2) |) | Т | Tr | nerm | no F | inn | nine | ette | (F3) | |
| ,001 | 10 2.4 | \vdash | <u></u> | 1 | | | | Τ | П | \dashv | Т | 7 | Τ | <u> </u> | T | Π | (, <u>_</u> / | | + | П | | T | | .p.,p. | Т | <u>, o,</u> | \Box |
| Mode | I / Range | 0.2-2ul | 1-10ul | 2-20ul | 2-20ul(µTip) | 5-50ul | 5-500l(µ lip) 10-100ul | 20-200ul | 30-300ul | 100-1000ul | 0.2-2ul(µTip) | 1-10ul(uTip) | 1 2 20ul | 2-20ul(µTip) | 2-50ul | 5-50ul(µTip) | 10-100ul | 20-2000l 100-1000ul | 0.2-2ul | 1-10nl | 1-10ul(µTip) | 2-20ul | 2-20ul(µTip) | 2-50ul | 5-50ul(µTip) | 10-100ul 20-200ul | 100-1000ul |
| Volume | Product No. | | | | | | | | | | Ĺ | Jnive | ersa | al Fit | Filt | er Ti | ips . | | | | | | | | | | |
| 0.5–10ul | 23-0010S | $\sqrt{}$ | $\sqrt{}$ | | | | | | | | 1 | | | | | | | | | $ \sqrt{} $ | | | | | | | |
| 0.5-1001 | 23-0011S | $ \sqrt{} $ | $\sqrt{}$ | | | | | | | | √ | | | | | | | | | $ \sqrt{} $ | | | | | | | |
| 1–20ul | 23–0020S 23–0021S | | | | | | | | | | | + | √ √ | | \vdash | | + | | + | | | $\frac{\sqrt{ }}{\sqrt{ }}$ | | | | + | |
| 1–50ul | 23-0050S 23-0051S | | | | | √ √ | √ √ | \[\sqrt{\sq}\}}}\sqrt{\sq}}\epsilon}\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sq}}\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sqrt{\sq}}}}}\sqrt{\sq}}}}}}\sqrt{\sqrt{\sq}}\sqrt{\sq}\sq}\sq}\sqrt{\sq}\sq}\sq\sintiq}\signt{\sq}\sq\sint{\sq}\signt{\sq}\signt{\sq}\sq}\signt{\sq}\signt{\sq}\signt{\sq}\signt{\sq}\signt{\sq}\signt{\sq}\si | | | + | | ļ | | √ √ | | $\sqrt{}$ | | H | | | | | √ √ | ۷ | | |
| 1–100ul | 23–0100S 23–0101S | | | | | √ | √ √ | \ | | | | | ļ | Ŧ | √ | | \[\langle \] | | F | | | | | √ | v | | |
| 1–300ul | 23-0200S | | | | | √ | √ | √ | √ _ | | | + | | ŧ | √ | | √ √ √ | | t | | | | | √ | ٧ | √ √ | |
| | 23-0201S | | | _ | \square | | + | + | √ | _ | 4 | + | - | + | \vdash | | √ | | _ | | | 4 | | Н | _ | | _ |
| 100- 1000ul | 23–1000S 23–1001S | | | + | \vdash | | + | + | Н | <u>√</u> | + | + | + | + | ╄ | | _ | √ | + | Н | | 4 | | \vdash | _ | + | √ |
| 100001 | 23-10013 | | | | | | | | Ш | √ | | | | | | | | 1 | | | | | | | | | ٧ . |
| Pipet | te Brand | | | BIC | HIT | Proli | ne | | | | BIO | HIT | Pro | oline | (P | lus) | | | | | | BR. | ΑN | D | | | |
| Mode | I / Range | 0.1-2.5ul | 0.5-10ul | 2-20ul | 5-50ul | 10-100ul | 20-200ul | 50-300ul | 100-1000ul | 0.1-3ul | 0.5-10ul | 2-20ul | 5-50ul | 10-100ul | 20-20011 | 30-30011 | 100-1000ul | 0.1-1ul | 0.1-2.5ul | 0.5-1001 | 2-2011 | | 5-50ul | 10-100ul | 20-200ul | 30-300ul | 100-1000ul |
| Volume | Product No. | ۳ | | " | ۵, | | '' | ۵, | | | | Jnive | | | | er Ti | | 10 | 10 | 1 | 7 (| ٠ ' | ۷, | | 1,, | (,) | |
| | 23–0010S | $\sqrt{}$ | | | | | | | | √ [| √ | ۷ | ., 00 | ا ۱ ا ا ا | | J. 11 | . ۲۵ | \int | | | - | | | | | | |
| 0.5–10ul | 23-0011S | \(| \(\frac{1}{}\) | | | | | | | <u>√</u> | $\sqrt{}$ | | | | | | | 1 | 1 | 1 | - | 1 | | | | | |
| 1 00 1 | 23-0020S | | | $\sqrt{}$ | | | | | | | | $\sqrt{}$ | | | | | | Ť | | | | 1 | | | | | |
| 1–20ul | 23-0021S | | | $\sqrt{}$ | | | | | | | | $\sqrt{}$ | | | | | | | | | | | | | | | |
| 1–50ul | 23-0050S | | | | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | | | | | $\sqrt{}$ | | | - √ | | | | | | ٦ | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | |
| 1-30ui | 23-0051S | | | | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | | | | | $\sqrt{}$ | _ √ | | - √ | | | | | | 1 | $\overline{}$ | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | |
| 1–100ul | 23-0100S | | | | | $\sqrt{}$ | <u>√</u> | $\sqrt{}$ | | | | | | | V | _ √ | | | | | | 1 | | $\sqrt{}$ | $\sqrt{}$ | $\sqrt{}$ | |
| | 23-01018 | | | | | √ | √ | √ | | | | | | √ | √ | - √ - _ | - | | | | | 1 | | √ | √ | √ | |
| 1-300ul | 23-0200S | | | | | | √ | √ | | | | | | | 1 | - √ - __ | - | | | | | + | | | √ | √ | |
| | 23–0201S 23–1000S | | | | | | √ | $\sqrt{}$ | _ | | | | | | √ | - √ | | + | | | | + | | | √ | $\sqrt{}$ | |
| 100- 1000ul | 23–1000S 23–1001S | | | | | | | | ν √ | | | | | | | | √ √ | + | | | | + | | | | | ν √ |
| 100001 | 20-10010 | | | | | | | | -4 | | | | | | | | ٧ | | | | | | | | | | ~ |

Rainin LTS Pipette Tips





>> Features

- Tips are made of medical grade polypropylene, and filters are made of highly hydrophobic **UHMWPE**
- 20µl, 200µl, and 1000µl tips
- Filters prevent cross-contamination and damage to samples, ensure accurate and precise sample aspiration, and protect the user by preventing sample liquid or vapor from entering the pipette body
- Low-retention tips undergo a coating-free treatment and will not contaminate samples.
- User-friendly labels on tip racks enable convenient usage
- DNase & RNase Free, non-pyrogen
- Sterile SAL 10⁻⁶



| Cat # | Description | Unit |
|----------|---|-----------------------------|
| R23-0020 | 20µl tip, rack, with filter, sterile | 96 Tips/Rack, 50 Racks/Case |
| R23-0021 | 20μl tip, rack, sterile, with filter, low retention | 96 Tips/Rack, 50 Racks/Case |
| R23-0200 | 200µl tip, rack, with filter, sterile | 96 Tips/Rack, 50 Racks/Case |
| R23-0201 | 200µl tip, rack, sterile, with filter, low retention | 96 Tips/Rack, 50 Racks/Case |
| R23-1000 | 1000µl tip, rack, with filter, sterile | 96 Tips/Rack, 50 Racks/Case |
| R23-1001 | 1000µl tip, rack, sterile, with filter, low retention | 96 Tips/Rack, 50 Racks/Case |

Non-filtered tips also available.

Automation Tips





- Made of medical grade medical grade polypropylene /black conductive medical grade polypropylene
- Compatible with Beckman, Tecan, Hamilton, and Agilent workstations
- Clear/black conductive tips, filtered/non-filtered tips, sterile/non-sterile tips low-retention/non low-retention tips available
- DNase & RNase Free, non-pyrogen
- Electron beam sterilization, SAL 10⁻⁶



| Cat # | Volume | Matched Workstations | Features | Length mm | Unit |
|-------------|--------|--|--|---------------------|---|
| AB23-0020S | 20µl | Beckman FX/NX | Clear, filtered, sterile, green racks | Standard (38.20) | 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case |
| AB23-0050S | 50µl | Beckman FX/NX | Clear, filtered, sterile, blue racks | Standard (51.20) | 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case |
| AB23-0250S | 250µl | Beckman FX/NX | Clear, filtered, sterile, purple racks | Standard (51.20) | 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case |
| AT23-0020S | 20µl | Tecan Genensis Freedom/ Freedom Evo | Clear, filtered, sterile | Standard (35.35) | 96 Tips/Rack, 24 Racks/Case |
| AT23-0020CS | 20μΙ | Tecan Genensis Freedom/ Freedom Evo | Conductive, filtered, sterile | Standard (35.35) | 96 Tips/Rack, 24 Racks/Case |
| AT23-0050S | 50µl | Tecan Genensis Freedom/ Freedom Evo | Clear, filtered, sterile | Standard (58) | 96 Tips/Rack, 24 Racks/Case |
| AT23-0050CS | 50µl | Tecan Genensis Freedom/ Freedom Evo | Conductive, filtered, sterile | Standard (58) | 96 Tips/Rack, 24 Racks/Case |
| AT23-0200S | 200μΙ | Tecan Genensis Freedom/ Freedom Evo | Clear, filtered, sterile | Standard (58.5) | 96 Tips/Rack, 24 Racks/Case |
| AT23-0200CS | 200μΙ | Tecan Genensis Freedom/ Freedom Evo | Conductive, filtered, sterile | Standard (58.5) | 96 Tips/Rack, 24 Racks/Case |
| AT23-1000S | 1000µl | Tecan Genensis Freedom/ Freedom Evo | Clear, filtered, sterile | Standard (95.5) | 96 Tips/Rack, 16 Racks/Case |
| AT23-1000CS | 1000µl | Tecan Genensis Freedom/ Freedom Evo | Conductive, filtered, sterile | Standard (95.5) | 96 Tips/Rack, 16 Racks/Case |
| AH23-0050S | 50µl | Hamilton Star/Starlet/Nimbus | Clear, filtered, sterile | Standard (50.4) | 96 Tips/Rack, 24 Racks/Case |
| AH23-0050CS | 50µl | Hamilton Star/Starlet/Nimbus | Conductive, filtered, sterile | Standard (50.4) | 96 Tips/Rack, 24 Racks/Case |
| AH23-0300S | 300µl | Hamilton Star/Starlet/Nimbus | Clear, filtered, sterile | Standard (60) | 96 Tips/Rack, 24 Racks/Case |
| AH23-0300CS | 300µl | Hamilton Star/Starlet/Nimbus | Conductive, filtered, sterile | Standard (60) | 96 Tips/Rack, 24 Racks/Case |
| AH23-1000S | 1000µl | Hamilton Star/Starlet/Nimbus | Clear, filtered, sterile | Standard (95) | 96 Tips/Rack, 16 Racks/Case |
| AH23-1000CS | 1000µl | Hamilton Star/Starlet/Nimbus | Conductive, filtered, sterile | Standard (95) | 96 Tips/Rack, 16 Racks/Case |
| AA21-0030S | 30µl | Agilent | Clear, sterilized | Standard (26.04) | 384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case |
| AA23-0250S | 250µl | Agilent | Clear, filtered, sterilized | Standard (55.50) | 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case |

Transfer Pipets







>> Features


- Polyethylene pipettes will not shatter or break
- An excellent alternative to Pasteur pipettes and rubber bulbs
- Evenly distributed thickness ensures uniform suction force
- Gamma radiation sterilized option available
- Available wrapped individually (sterile) or in bulk package (non-sterile)

Molded graduations

Transparent wall permits easy viewing of pipette contents







Dimensions = mm


| Cat # | Length | Graduated To | Sterile | Package | Unit | Total |
|-----------|---------|--------------|---------|----------------------|-------------------------|------------|
| 30-0135 | 162.0mm | 1ml | No | Bulk | 500 /Pack, 4 Packs/Case | 2000 /Case |
| 30-0135A1 | 162.0mm | 1ml | Yes | Individually Wrapped | 500 /Pack, 4 Packs/Case | 2000 /Case |
| 30-0138 | 162.0mm | 3ml | No | Bulk | 500 /Pack, 4 Packs/Case | 2000 /Case |
| 30-0138A1 | 162.0mm | 3ml | Yes | Individually Wrapped | 500 /Pack, 4 Packs/Case | 2000 /Case |
| 30-0238A1 | 183.0mm | 3ml | Yes | Individually Wrapped | 500 /Pack, 4 Packs/Case | 2000 /Case |


Solution Basins





- Made of PVC, polystyrene (PS), or strong polypropylene (PP)
- PVC basins have transparent walls to permit easy viewing of contents
- Polypropylene basins are chemically resistant to alcohols and mild organic solvents
- PS one-channel basins are available in three sizes: 25ml, 55ml, and 100ml; 25ml and 55ml basins have a V-shaped bottom design
- PVC one-channel basins have a capacity of 55ml and a V-shaped bottom design
- PP 12-channel basins are double sided: the 12-channel side has a capacity of 5ml per channel; the other side (single-channel) has a capacity of 50ml
- PP 12-channel basins are available in assorted colors: clear, blue, green, pink, and yellow
- PS basins are individually packaged (sterile); PVC and PP basins are bulk packaged (non-sterile)





Dimensions = mm


| Cat # | Volume | Material | Color | Sterile | Package | Unit | Total |
|---------|--------|----------|-----------------|---------|----------------------|-----------------------|-----------|
| 25-1025 | 25ml | PS | White | Yes | Individually Wrapped | 1 /Bag, 100 Bags/Case | 100 /Case |
| 25-1202 | 50ml | PP | Assorted Colors | No | Bulk | 25 /Bag, 4 Bags/Case | 100 /Case |
| 25-0052 | 55ml | PVC | Clear | No | Bulk | 100 /Bag, 8 Bags/Case | 800 /Case |
| 25-0051 | 55ml | PS | White | Yes | Individually Wrapped | 1 /Bag, 100 Bags/Case | 100 /Case |
| 25-1100 | 100ml | PS | White | Yes | Individually Wrapped | 1 /Bag, 50 Bags/Case | 50 /Case |

Beakers






- Made of strong, medical grade polypropylene, autoclavable
- Tri-corner design provides three dripless pouring spouts
- Six different sizes meet different needs



| Cat # | Volume | Sterile | Unit | Total |
|---------|---------|---------|----------------------------------|-----------|
| 93-0050 | 50ml | No | 25 Pieces/Sleeve, 4 Sleeves/Case | 100 /Case |
| 93-0100 | 100ml | No | 25 Pieces/Sleeve, 4 Sleeves/Case | 100 /Case |
| 93-0250 | 250ml | No | 25 Pieces/Sleeve, 4 Sleeves/Case | 100 /Case |
| 93-0400 | 400ml | No | 25 Pieces/Sleeve, 4 Sleeves/Case | 100 /Case |
| 93-0800 | 800ml | No | 25 Pieces/Sleeve, 4 Sleeves/Case | 100 /Case |
| 93-1000 | 1,000ml | No | 25 Pieces/Sleeve, 4 Sleeves/Case | 100 /Case |



 ${\sf Dimensions} = {\sf mm}$

Reagent Bottles





- Made of durable medical grade polypropylene or high-density polyethylene
- Uniform wall thickness and strength, impact-resistant, great chemical resistance to most acid, alkali, and alcohol
- Polypropylene bottles are autoclavable
- Wide-mouth design enables the superior sealing performance, and the convenience for filling liquid
- Lineless caps provide the ultimate in leakproof protection without the use of liner that can wrinkle, cause leaks, or contaminate your reagents
- Amber/browm bottles reduce UV light transmissions to protect light-sensitive liquids; Clear/natural bottles are convenient to observe inner items
- Various volumes are available: 8ml, 15ml, 30ml, 60ml, 125ml, 250ml, 500ml, and 1000ml
- DNase & RNase free, non-pyrogenic
- Sterile bottles also available



| Cat # | Color | Cat # | Color | Volume | Material | Sterile | Unit |
|----------|---------|----------|-------|--------|-------------------------|---------|------------------------|
| 04-0008U | Clear | 04-1008U | Amber | 8ml | PP | No | 100 /Bag, 10 Bags/Case |
| 04-0015U | Clear | 04-1015U | Amber | 15ml | PP | No | 100 /Bag, 10 Bags/Case |
| 04-0030U | Clear | 04-1030U | Amber | 30ml | PP | No | 100 /Bag, 10 Bags/Case |
| 04-0060U | Clear | 04-1060U | Amber | 60ml | PP | No | 50 /Bag, 10 Bags/Case |
| 04-0125U | Clear | 04-1125U | Amber | 125ml | PP | No | 25 /Bag, 10 Bags/Case |
| 04-0250U | Clear | 04-1250U | Amber | 250ml | PP | No | 20 /Bag, 10 Bags/Case |
| 04-0500U | Clear | 04-1500U | Amber | 500ml | PP | No | 10 /Bag, 10 Bags/Case |
| 04-0100U | Clear | 04-1100U | Amber | 1000ml | PP | No | 5 /Bag, 20 Bags/Case |
| 04-2008 | Natural | 04-3008 | Brown | 8ml | Bottles: HDPE, caps: PP | No | 1000 /Case |
| 04-2015 | Natural | 04-3015 | Brown | 15ml | Bottles: HDPE, caps: PP | No | 1000 /Case |
| 04-2030 | Natural | 04-3030 | Brown | 30ml | Bottles: HDPE, caps: PP | No | 400 /Case |
| 04-2060 | Natural | 04-3060 | Brown | 60ml | Bottles: HDPE, caps: PP | No | 650 /Case |
| 04-2125 | Natural | 04-3125 | Brown | 125ml | Bottles: HDPE, caps: PP | No | 350 /Case |
| 04-2250 | Natural | 04-3250 | Brown | 250ml | Bottles: HDPE, caps: PP | No | 200 /Case |
| 04-2500 | Natural | 04-3500 | Brown | 500ml | Bottles: HDPE, caps: PP | No | 100 /Case |
| 04-2100 | Natural | 04-3100 | Brown | 1000ml | Bottles: HDPE, caps: PP | No | 50 /Case |


Test Tubes





- Tubes: made of medical grade polystyrene or polypropylene
- Caps: made of medical grade polyethylene
- Tubes are available in two sizes: 12×75 mm and 13×100 mm





| Cat # | Description | Size(mm) | Material | Color | Volume | Sterile | Unit | Total |
|---------|-------------|----------|----------|-------|--------|---------|-------------------------|------------|
| 12-1275 | Tube | 12×75 | PS | Clear | 5ml | No | 500 /Bag, 4 Bags/Case | 2000 /Case |
| 12-2275 | Tube | 12×75 | PP | Clear | 5ml | No | 500 /Bag, 4 Bags/Case | 2000 /Case |
| 12-1300 | Tube | 13×100 | PS | Clear | 8ml | No | 250 /Bag, 8 Bags/Case | 2000 /Case |
| 12-1299 | Cap | 12x13.5 | PE | Blue | - | No | 1,000 /Bag, 2 Bags/Case | 2000 /Case |
| 12-1399 | Cap | 13x17.5 | PE | Red | - | No | 1,000 /Bag, 2 Bags/Case | 2000 /Case |

Slide Storage Boxes






>> Features

Microscope Slides

- Made of high quality, highly transparent glass
- Slide surface highly polished
 Ground edges prevent chipping and personal injury
- Slides with frosted area are convenient for labeling and storage Slides without frosted area have a larger usable surface area

Slide Storage Boxes

- Made of strong ABS material
- Available in 25-place or 100-place sizes
- Traditional cork or foam lining available
- 5 assorted colors available: red, yellow, green, blue, and white
- Store 25mm × 75mm microscope slides
- Each slot is numbered to correspond to the slide inventory sheet inside the lid
- 100-place slide boxes are secured with a rust-resistant clasp and hinge pin



Dimensions = mm

Microscope Slides

| Cat # | Material | Туре | Size (mm) | Sterile | Unit | Total |
|-------|-------------|--------------------|-------------|---------|-------------------------|------------|
| 7101 | Clear Glass | Plain Slides | 25.4 x 76.2 | No | 72 /Pack, 50 Packs/Case | 3600 /Case |
| 7105 | Clear Glass | Single Frosted End | 25.4 x 76.2 | No | 72 /Pack, 50 Packs/Case | 3600 /Case |

Slide Storage Boxes


| Cat # | Place | Lining | Color | Sterile | Unit | Total | |
|---------|-----------|--------|----------|---------|---------------------|----------|--|
| | | | | | | | |
| 41-5025 | 25-place | Foam | Assorted | No | 5 /Bag, 4 Bags/Case | 20 /Case | |
| 41-6025 | 25-place | Cork | Assorted | No | 5 /Bag, 4 Bags/Case | 20 /Case | |
| 41-5100 | 100-place | Foam | Assorted | No | 5 /Bag, 6 Bags/Case | 30 /Case | |
| 41-6100 | 100-place | Cork | Assorted | No | 5 /Bag 6 Bags/Case | 30 /Case | |

Deep Well Plates






- Made of medical grade polypropylene, autoclavable
- Chemically resistant to alcohols and mild organic solvents
- 0.5ml, 1.0ml, 1.6ml, 2.2ml and 10ml volumes available
- Square and round wells, V and U bottom available
- DNase & RNase free, pyrogen free
- Can be sealed with sealing films or silicone sealing mats
- Confirms with SBS standards for automation
- Compatible with automated nucleic acid extraction equipment
- Sterilized, individual wrapped, and with lid type optional



| Cat # | Volume | Description | Сар | Sterile | Unit |
|------------|--------|--|-----|---------|--------------------------------|
| 02-8005 | 0.5ml | Square 96 wells, V bottom, KingFisher Flex system | No | No | 10 Pieces/Pack,100 Pieces/Case |
| 02-2010-50 | 1.0ml | Round 96 wells, U bottom, | No | No | 10 Pieces/Pack, 50 Pieces/Case |
| 02-1016 | 1.6ml | Square 96 wells, U bottom | No | No | 24 Pieces/Pack, 96 Pieces/Case |
| 02-1022 | 2.2ml | Square 96 wells, U bottom | No | No | 24 Pieces/Pack, 96 Pieces/Case |
| 02-6022 | 2.2ml | Square 96 wells, U bottom, profile concave | No | No | 24 Pieces/Pack, 96 Pieces/Case |
| 02-7022 | 2.2ml | Square 96 wells, V bottom, KingFisher Flex system | No | No | 10 Pieces/Pack, 50 Pieces/Case |
| 02-9010 | 10ml | Square 24 wells, V bottom, KingFisher Flex system | No | No | 10 Pieces/Pack, 50 Pieces/Case |
| 02-0096 | / | 96-well magnet set | / | No | 2 Pieces/Pack, 50 Pieces/Case |
| 02-0008 | / | 8-strip magnet set | / | No | 2 Pieces/Pack, 100 Pieces/Case |
| 02-1096 | / | 96 square silicone sealing mats, for 96 square deep well plates | / | No | 10 Pieces/Pack, 50 Pieces/Case |
| 02-1196 | / | 96 square silicone sealing mats, can be punctured, for 96 square deep well plates | / | No | 10 Pieces/Pack, 50 Pieces/Case |
| 02-2196 | / | 96 round silicone sealing mats, can be punctured, for 1.0ml plates | / | No | 10 Pieces/Pack, 50 Pieces/Case |
| 02-2296 | / | 96 round silicone sealing mats, can be punctured, for 2.0ml plates | / | No | 10 Pieces/Pack, 50 Pieces/Case |
| 02-3024 | / | 24 square silicone sealing mats, for 24 square deep well plates 10ml | / | No | 10 Pieces/Pack, 50 Pieces/Case |

^{*} Sterilized, individual wrapped, and with lid type optional.

Sealing Films





- Made of High quality polypropylene, transparent polyester, hard polyester, impermeable aluminum foil, vinyl
- Well with not only PCR, qPCR, ELISA, inoculation and fluorescence experiment but also general and high throughput operation
- Applied to ELISA, inoculation, sample storage and microorganism cultivation sterile/ non-sterile
- Multiple temperature ranges are available, -80 $^{\circ}$, -70 $^{\circ}$, -40 $^{\circ}$, -20 $^{\circ}$, +80 $^{\circ}$, 90 $^{\circ}$, +100 $^{\circ}$, +120 $^{\circ}$
- Colored sealers with five options are available, green/purple/orange/red/blue
- Film-sealing paddle and plate roller are optional



| Cat# | Material | Temperature range | Application | Sterile | Unit | Total |
|---------|--------------------------------|-------------------|--|---------|---|-----------|
| 61-0010 | Polypropylene | -40℃ ~ +120℃ | ELISA, inoculation | No | 100 Pieces/Box | 100 /Case |
| 61-0020 | Impermeable aluminum foil | -80℃ ~+120℃ | PCR, cold storage of sample | No | 100 Pieces/Box | 100 /Case |
| 61-0040 | Hard polyester | -40℃ ~ +120℃ | PCR, qPCR and sample storage | No | 100 Pieces/Box | 100 /Case |
| 61-0041 | Hard polyester | -40℃ ~+120℃ | PCR, qPCR and sample storage | No | 100 pieces/Box | 100 /Case |
| 61-0051 | Thick porous hydrophobic rayon | -20℃ ~+80℃ | biological culture | Yes | 50 Pieces/Box | 50 /Case |
| 61-0060 | Thick vinyl | -40℃ ~ +90℃ | automatic sampling | No | 100 Pieces/Box | 100 /Case |
| 61-0080 | DMSO resistant sealing foil | -80℃ ~+120℃ | PCR, photosensitive test and mechanical operation | No | 100 Pieces/Box | 100 /Case |
| 61-0120 | Hard polyester | -40℃ ~ +120℃ | ELISA, inoculation | Yes | 2 Rolls/Box, 100 Pieces/Roll, equip with dispenser | 200 /Case |
| 61-0152 | Thick porous hydrophobic rayon | -20℃ ~+80℃ | biological culture | Yes | 2 Rolls/Box, 50 Pieces/Roll, equip with dispenser | 100 /Case |
| 61-0210 | Polyolefin | -70℃ ~ +100℃ | qPCR, sample storage and protein crystallization | No | 100 Pieces/Box | 100 /Case |
| 61-0303 | Transparent polyester | -40℃ ~+120℃ | PCR, qPCR and sample storage | No | 400 Pieces/Box | 400 /Case |
| 61-0304 | Aluminum foil | -80℃ ~ +120℃ | PCR, sample storage photosensitive test, automatic operation | No | 300 Strip/Box | 300 /Case |






- Available in powder free latex and powder free nitrile
- Textured fingertips for an improved grip Beaded cuff design helps users don gloves efficiently
- Non-sterile






| Cat # | Material | Size | Color | Unit | Total |
|---------|----------|------|-------------------|--------------------------|------------|
| 97-1111 | Latex | XS | Natural Off-White | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-1112 | Latex | S | Natural Off-White | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-1113 | Latex | М | Natural Off-White | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-1114 | Latex | L | Natural Off-White | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-1115 | Latex | XL | Natural Off-White | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-6111 | Nitrile | XS | Blue | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-6112 | Nitrile | S | Blue | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-6113 | Nitrile | М | Blue | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-6114 | Nitrile | L | Blue | 100 /Pack, 10 Packs/Case | 1000 /Case |
| 97-6115 | Nitrile | XL | Blue | 100 /Pack, 10 Packs/Case | 1000 /Case |



Screw Cap Microtubes



- Made of medical grade polypropylene
- Can be repeatedly frozen and thawed
- Universal screw threads for use
- Tubes fit most common rotors
- Tubes fit standard 1-inch and 2-inch
- Autoclavable to 121°C and freezable to -86°C



| Cat# | Туре | Volume | Height | Color | Bottom | Sterile | Unit | Total |
|----------|-------|--------|--------|----------|---------------|---------|-----------------------------|------------|
| 81-0053 | Tubes | 0.5ml | 30.0 | Clear | Conical | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0054 | Tubes | 0.5ml | 31.0 | Clear | Self-Standing | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0055 | Tubes | 0.5ml | 47.5 | Brown | Self-Standing | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0056 | Tubes | 0.5ml | 47.5 | Clear | Self-Standing | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0153 | Tubes | 1.5ml | 46.5 | Clear | Conical | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0154 | Tubes | 1.5ml | 47.5 | Clear | Self-Standing | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0203 | Tubes | 2.0ml | 46.1 | Clear | Conical | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0204 | Tubes | 2.0ml | 47.5 | Clear | Self-Standing | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0205 | Tubes | 2.0ml | 47.5 | Brown | Self-Standing | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0000 | Caps | - | 9.7 | Clear | - | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0001* | Caps | - | 9.7 | Pink | - | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0002* | Caps | - | 9.7 | Green | - | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0003* | Caps | - | 9.7 | Blue | - | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0004* | Caps | - | 9.7 | Yellow | - | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0005* | Caps | - | 9.7 | Assorted | - | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |
| 81-0006 | Caps | - | 9.7 | Brown | - | No | 500 /Bag/Pack,10 Packs/Case | 5000 /Case |


^{*} Call for details.

Cryogenic Vials





- Made of medical grade polypropylene
- Autoclavable
- Can be repeatedly frozen and thawed
- Universal screw threads
- Vials fit standard 1-inch and 2-inch freezer boxes
- Gamma radiation sterilized
- DNase & RNase free, endotoxin free, and human DNA free
- 81-8204 2ml self-standing vials have a white writing area and black graduations



| Cat # | Volume | Color | Height | Bottom | Writing Area | Sterile | Unit | Total |
|----------|--------|-------|--------|---------------|--------------|---------|-------------------------------------|------------|
| 81-7053* | 0.5ml | Clear | 30.0 | Conical | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7054 | 0.5ml | Clear | 31.0 | Self-Standing | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7055 | 0.5ml | Brown | 47.5 | Self-Standing | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7056 | 0.5ml | Clear | 47.5 | Self-Standing | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7153 | 1.5ml | Clear | 46.5 | Conical | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7154 | 1.5ml | Clear | 47.5 | Self-Standing | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7203 | 2.0ml | Clear | 46.1 | Conical | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7204 | 2.0ml | Clear | 47.5 | Self-Standing | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-7205 | 2.0ml | Brown | 47.5 | Self-Standing | No | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |
| 81-8204 | 2.0ml | Clear | 47.5 | Self-Standing | Yes | Yes | 50 /Bag, 10 Bags/Pack, 4 Packs/Case | 2000 /Case |

^{*} Call for details.

Traditional Cryogenic Vials



>> Features

- Made of medical grade polypropylene; Autoclavable; and can be repeatedly frozen and thawed
- 5ml vials are available with external screw threads, and 2ml vials are available with both external and internal screw threads. External thread screw vials' caps have no O-ring, which reduces the possibility of contamination
- Fit 9×9 (81-well) standard cryogenic boxes
- Tubes with black graduation printing and white writing area
- Gamma radiations sterilized, RNase and DNase free, endotoxin free and foreign DNA free

2ml External

Thread

• Pre-set 1D barcodes and numbers on sides for immediate scanning

Working temperature: stable from -196°C to 121°C Suitable for liquid nitrogen freezing(vapor phase)

88-3202S

Dimensions=mm

2ml External

Thread-S Series

5ml External

Thread

| Cat # | Type | Thread Type | Volume | Bottom | Side Barcode | Writing Area | Sterile | Unit | Total |
|----------|-----------|-------------|----------|---------------|--------------|--------------|---------|----------------------------|----------------|
| 88-321X | Tubes | Internal | 2ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-320X | Tubes | External | 2ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-350X | Tubes | External | 5ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-621X | Tubes | Internal | 2ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-620X | Tubes | External | 2ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-650X | Tubes | External | 5ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-600X | Markers | - | - | - | - | - | - | 100 Sets/Bag, 30 Bags/Case | 3000 Sets/Case |
| 88-320XS | Tubes | External | 2ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-350XS | Tubes | External | 5ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-620XS | Tubes | External | 2ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-650XS | Tubes | External | 5ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| X—Col | or of Can | s:0-White | 1-Red 2- | -Green 3-Blu | e 4-Yellow 5 | 5-Purple, 7- | Clear | 80-6006-Assorted Co | lor |

2ml Internal

Thread

5ml External

Thread-S Series

CryoKING Cryogenic Vials






>> Features

- Made of medical grade polypropylene, can be repeatedly frozen and thawed
- Vials have a white writing area and black graduations; 6 different colored caps are available
- Vials with pre-set Datamatrix 2D barcodes are available, and the 2D barcodes can be customized
- Use the microscopic imaging technique to pre-set 2D barcodes
- Gamma radiation sterilized, DNase & RNase free, endotoxin free, foreign DNA free
- Pre-set 1D barcodes and human readable numbers on sides for visual identification, and can be read by various scanners

Specifications

- Working temperature: stable from -196°C to 121°C
- Suitable for liquid nitrogen freezing(vapor phase)



Dimensions = mm

| Cat # | Cat# (Upgraded, with Side Barcode) | Thread Type | Volume | Bottom | Bottom Barcode | Writing Area | Sterile | Unit | Total |
|---------|--|-------------|--------|---------------|-------------------|-----------------|---------|---------------------------|----------------|
| 88-005X | 88-605X | External | 0.5ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-010X | 88-610X | External | 1.0ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-015X | 88-615X | External | 1.5ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-020X | 88-720X | External | 2.0ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-021X | 88-721X | Internal | 2.0ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| 88-050X | 88-750X | External | 5.0ml | Self-Standing | No | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| - | 88-905X | External | 0.5ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| _ | 88-910X | External | 1.0ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| - | 88-915X | External | 1.5ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| _ | 88-920X | External | 2.0ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| - | 88-921X | Internal | 2.0ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |
| _ | 88-950X | External | 5.0ml | Self-Standing | Yes | Yes | Yes | 25 Sets/Bag, 40 Bags/Case | 1000 Sets/Case |

CryoKING Combo





>> Features

- All you need in one package:
 - CryoKING Vials with pre-set 2D barcodes
 - Pre-assembled inside the CryoKING Boxes
 - Pre-sterilized
- Certified: foreign DNA free, DNase & RNase free, endotoxin free
- Rapidly scans and identifies an entire box full of vials
- Pre-set 1D barcodes and numbers on side for immediate scanning or visual identification.

Specifications

- Working temperature: stable from -196°C to 121°C
- Suitable for liquid nitrogen freezing(vapor phase)


Boxes: 9*9 or 10*10 wells for bottom coded vials



Vials: Multiple codes for efficient sample management

| Cat # | Volume | Bottom barcode | Side Code | Wells | Thread Type | Dimensions(mm) | | Unit | Total | |
|---|-------------------------|-------------------|-----------|----------|-------------|----------------|--------------|--------------|--------------|--|
| 89-905X | 0.5ml Vials + 1 inBox | Yes | Yes | 100-well | External | 133x133x36 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| 89-910X | 1.0ml Vials + 2 inBox | Yes | Yes | 100-well | External | 133x133x52 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| 89-915X | 1.5ml Vials + 2 inBox | Yes | Yes | 100-well | External | 133x133x52 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| 89-922X | 2.0ml Vials + 2 inBox | Yes | Yes | 100-well | External | 136x136x52 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| 89-923X | 2.0ml Vials + 2 inBox | Yes | Yes | 100-well | Internal | 136x136x52 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| 89-920X | 2.0ml Vials + 2 inBox | Yes | Yes | 81-well | External | 133x133x52 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| 89-921X | 2.0ml Vials + 2 inBox | Yes | Yes | 81-well | Internal | 133x133x52 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| 89-950X | 5.0ml Vials + 3.75 inBo | x Yes | Yes | 81-well | External | 133x133x95 | 6 Bags/Pack, | 2 Packs/Case | 12 Bags/Case | |
| XCap Color : 0-White,1-Red,2-Green,3-Blue,4-Yellow,5-Purple,7-Clear | | | | | | | | | | |

Automated Product Series of SBS Standards-Internal Thread





>> Features

- 2D barcodes on 96-well format cryogenic vials are easy to be read
- 2D barcode and human-readable numbers on the side of SBS format racks; 2D barcode at the bottom of racks, cryogenic vials and racks can be traced at the same time
- Simple and powerful mechanical locking function, slide the buttons on both sides to unlock
- Laser etched 2D barcodes, resistant to DMSO and other organic solvents
- SBS format rack fits 96-well format cryogenic vials perfectly for automatic operations

Specifications

- Working temperature: stable from -196°C to 121°C
- Suitable for liquid nitrogen freezing(vapor phase)





96-well SBS Cryogenic Vials

| Cat # | Thread Type | Volume | Bottom | Bottom Barcode | Writing Area | Sterile | Unit | Total |
|---------|--------------|-----------|-----------------|-------------------|-----------------|---------|---|----------------|
| 88-403X | Internal | 300µl | Self-Standing | Yes | No | Yes | 96 Sets/Bag, 10 Bags/Pack, 2 Packs/Case | 1920 Sets/Case |
| 88-405X | Internal | 500µl | Self-Standing | Yes | No | Yes | 96 Sets/Bag, 10 Bags/Pack, 2 Packs/Case | 1920 Sets/Case |
| 88-408X | Internal | 750µl | Self-Standing | Yes | No | Yes | 96 Sets/Bag, 10 Bags/Pack, 2 Packs/Case | 1920 Sets/Case |
| 88-410X | Internal | 1000µl | Self-Standing | Yes | No | Yes | 96 Sets/Bag, 10 Bags/Pack, 2 Packs/Case | 1920 Sets/Case |
| XCap (| Color : 0-Wh | ite.1-Red | l.2-Green.3-Blu | e.4-Yellov | v.5-Purpl | e | | |

SBS Format Cryogenic Racks

| Cat # | Vial Volume | Wells | Color | Dimensions(mm) | Unit | Total |
|---------|----------------|-----------------|-------|----------------|-------------------------|----------|
| 98-4127 | 300µl & 500µl | 96-well, square | White | 128x85x35 | 10 /Pack , 2 Packs/Case | 20 /Case |
| 98-4227 | 750µl & 1000µl | 96-well, square | White | 128x85x54 | 10 /Pack , 2 Packs/Case | 20 /Case |

SBS Format Combo

| Cat # | Description | Bottom Barcode | Side Barcode | Wells | Dimensions(mm) | Unit |
|-----------|-------------------|-------------------|--------------|-----------------|----------------|---|
| 89-423X | 300µl Vials+Rack | Yes | No | 96-well, square | 128x85x35 | 1 Set/Bag , 10 Bags/Pack , 2 Packs/Case |
| 89-425X | 500µl Vials+Rack | Yes | No | 96-well, square | 128x85x35 | 1 Set/Bag , 10 Bags/Pack , 2 Packs/Case |
| 89-428X | 750µl Vials+Rack | Yes | No | 96-well, square | 128x85x54 | 1 Set/Bag , 10 Bags/Pack , 2 Packs/Case |
| 89-421X 1 | 1000µl Vials+Rack | Yes | No | 96-well, square | 128x85x54 | 1 Set/Bag , 10 Bags/Pack , 2 Packs/Case |

X--Cap Color: 0-White, 1-Red, 2-Green, 3-Blue, 4-Yellow, 5-Purpl

Automated Products Series of SBS Standards-External Thread





>> Features


- Vials: Made of medical grade polypropylene, 2-Color injection molding technique. Four codes on one vial: Bottom: Laser etched 2D code and human-readable codes; Side: Laser etched 1D barcode and human-readable codes.
- Racks: Side: Laser etched 1D barcode, 2D code, and human-readable codes; Bottom: 2D code. Pre-set alphanumeric marks for sample identification. Can be used with CKS-1301 CryoKING scanner to manager samples.
- 100,000 grade clean room enables safe vials.
- Gamma radiations sterilized; DNase & RNase free, endotoxin free, foreign DNA free.
- Stable from -196°C to 121°C, can be stored in vapor phase liquid nitrogen.
- Compatible with Micronic cappers & decappers and sample management system.



| Cat # | Description | Unit |
|---------|--|---|
| 88-514X | CryoKING SBS vials, 1.4ml, external thread; Self-Standing; 2-color injection molding technique; Laser-etched 4 codes. Sterile. | 96 Sets/Bag, 10 Bags/Pack, 2 Packs/Cas |
| 98-5227 | CryoKING SBS format racks, available for 1.4ml external vials; 96-well; Square holes, White; Dimensions(mm): 128x85x56. Non sterile. | 10 /Pack, 2 Packs/Case |
| 89-514X | CryoKING SBS Combo, 1.4ml Vials + Rack. Dimensions(mm): 128x85x56. Sterile. | 1 Set/Bag, 10 Bags/Pack, 2 Packs/Case |

X—Cap Color: 1-Sky blue, 2-Lime green, 3-Biologix blue, 5-Lavender purple

Premium Cardboard Cryogenic Boxes





>> Features

- Made of durable cardboard with superior white outer coating
- Chemically resistant to alcohols and mild organic solvents
- Five assorted colors available: white, yellow, pink, blue, and green
- White Premium Cardboard Freezer Boxes have a glossy coating, other colors have a matte coating
- Water resistant coating for improved durability
- Superior white 90-1281 & 90-1200 cardboard cryogenic boxes without outer coating are also available
- Holds 0.5ml, 1.5ml, 2.0ml, 5ml and 10ml microtubes and 15ml or 50ml Centrifuge Tubes

- Working temperature: stable from -196°C to 121°C
- Suitable for liquid nitrogen freezing



| Cat # | Height | Wells | Description | Dimensions (mm) | Unit | Total |
|---------|--------|----------|--------------------------------------|-----------------|----------------------|-----------|
| | | | | | | |
| 90-1200 | 2 in | 100-well | Superior white without outer coating | 132 x 132 x 53 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-1281 | 2 in | 81-well | Superior white without outer coating | 132 x 132 x 53 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-1300 | 3 in | 100-well | Superior white without outer coating | 133 x 133 x 76 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-1381 | 3 in | 81-well | Superior white without outer coating | 133 x 133 x 76 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-2281 | 2 in | 81-well | White coating, plastic coating | 133 x 133 x 52 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-2200 | 2 in | 100-well | White coating, plastic coating | 133 x 133 x 52 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-5281 | 2 in | 81-well | Rainbow outer coating, matte coating | 133 x 133 x 52 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-5200 | 2 in | 100-well | Rainbow outer coating, matte coating | 133 x 133 x 52 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-2381 | 3 in | 81-well | White coating, plastic coating | 133 x 133 x 76 | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-2300 | 3 in | 100-well | White coating, plastic coating | 133 x 133 x 76 | 5 /Bag, 20 Bags/Case | 100 /Case |

ID-Color™ Cardboard Cryogenic Boxes







>> Features

- Made of durable cardboard with a waterproof coating and polypropylene dividers
- Connected lid design ensures that lids are never misplaced or mixed up
- Alphanumeric labeling featured on the lid, bottom, and four sides enables quick and easy identification of samples
- Holds 0.5ml/1.5ml/2.0ml/5.0ml microtubes and cryogenic vials

Specifications

• Suitable for liquid nitrogen freezing (after adding vent holes, and using cardboard dividers)






Straight grid divider design, resistant to low temperature

| Cat # | Height | Well | Divider Material | Dimensions (mm) | Unit | Total |
|---------|--------|----------|------------------|-----------------|---------------------|----------|
| 90-8100 | 1in | 100-well | PP | 133 x 133 x 36 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-8281 | 2in | 81-well | PP | 133 x 133 x 52 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-8200 | 2in | 100-well | PP | 133 x 133 x 52 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-8381 | 3.75in | 81-well | PP | 133 x 133 x 95 | 5 /Bag, 4 Bags/Case | 20 /Case |

PC Cryogenic Boxes





>> Features

- Made of tough polycarbonate, chemically resistant to alcohols and mild organic solvents
- Can be repeatedly frozen and thawed
- 5x5 (25-well), 9x9 (81-well), or 10x10 (100-well) boxes available
- Compatible with 0.5ml, 1.5ml, 2.0ml or 5.0ml microtubes and cryogenic vials
- Numbered grids for sample identification
- Drain holes and air vents minimize condensation
- Five assorted colors available

- Suitable for storage in freezers or liquid nitrogen
- Temperature range: stable from -196°C to 121°C










Numbers for convenient sample organization

98-series, bottom holes for bottom coded vials

| Cat # | Height | Wells | Color | Dimensions (mm) | Unit | Total |
|---------|---------|----------|-----------------|-----------------|---------------------|----------|
| 90-9009 | 2 in | 81-well | Assorted Colors | 133x133x52 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-9210 | 2 in | 100-well | Assorted Colors | 133x133x52 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-9250 | 2 in | 25-well | Assorted Colors | 75x75x52 | 5 /Bag, 6 Bags/Case | 30 /Case |
| 90-9381 | 3.75 in | 81-well | Assorted Colors | 133x133x95 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 98-0113 | 1 in | 100-well | Blue | 133x133x36 | 3 /Bag, 4 Bags/Case | 12 /Case |
| 98-0200 | 2 in | 100-well | Blue | 136x136x52 | 3 /Bag, 4 Bags/Case | 12 /Case |
| 98-0213 | 2 in | 100-well | Blue | 133x133x52 | 3 /Bag, 4 Bags/Case | 12 /Case |
| 98-0313 | 3.75 in | 81-well | Blue | 133x133x95 | 3 /Bag, 4 Bags/Case | 12 /Case |
| 98-1213 | 2 in | 81-well | Blue | 133x133x52 | 3 /Bag, 4 Bags/Case | 12 /Case |

PP Cryogenic Boxes






>> Features

- Made of polypropylene, chemically resistant to alcohols and mild organic solvents
- Can be repeatedly frozen and thawed
- 5x5 (25-well), 9x9 (81-well), or 10x10 (100-well) boxes available
- Compatible with 0.5ml, 1.5ml or 2.0ml microtubes and cryogenic vials
- Numbered grids for sample identification
- Drain holes and air vents minimize condensation
- Five assorted colors available

- Suitable for storage in freezers
- Temperature range: stable from -86°C to 121°C





Numbers for convenient sample organization

| Cat # | Height | Wells | Color | Dimensions (mm) | Unit | Total |
|---------|--------|----------|-----------------|-----------------|---------------------|----------|
| 90-9025 | 2 in | 25-well | Assorted Colors | 75 x 75 x 52 | 5 /Bag, 6 Bags/Case | 30 /Case |
| 90-9081 | 2 in | 81-well | Assorted Colors | 133 x 133 x 52 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-9100 | 2 in | 100-well | Assorted Colors | 141 x 154 x 55 | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-9200 | 2 in | 100-well | Assorted Colors | 133 x 133 x 52 | 5 /Bag, 4 Bags/Case | 20 /Case |

Freezer Racks



Sliding Drawer Freezer Racks

- Made of corrosion-resistant stainless steel or aluminium alloy
- Ideal for storage of 1 inch and 2 inch standard size cryogenic boxes
- Fits most standard upright freezers
- Sliding Drawer Freezer Racks include folded rails to reduce frost
- Flanging design reduces scratching
- Handle design enables easy removal from the shelf
- Colored labels for drawer/sample identification
- Can be customized

| Cat # | Material | Description | Capacity (D x H Configuration) | Dimensions(mm) (D x W x H) | Storage box size (mm) |
|---------|-------------------|------------------------------|-----------------------------------|-------------------------------|--------------------------|
| | | | | | |
| 99-7216 | Stainless steel | Sliding Drawer Freezer Racks | 4 x 4 | 564 x 141 x 229 | 133 x 133 x 52 |
| 99-7220 | Stainless steel | Sliding Drawer Freezer Racks | 4 x 5 | 564 x 141 x 285 | 133 x 133 x 52 |
| 99-1220 | Otali liess steel | Oliding Drawer Treezer Nacks | 4 X 3 | 304 X 141 X 263 | 100 X 100 X 02 |
| 99-7224 | Stainless steel | Sliding Drawer Freezer Racks | 4 x 6 | 564 x 141 x 341 | 133 x 133 x 52 |
| 96-7216 | Aluminium alloy | Sliding Drawer Freezer Racks | 4 x 4 | 564 x 141 x 231 | 133 x 133 x 52 |
| 00 7000 | Aluminium allay | Sliding Drawer Freezer Racks | 4 5 | 504 444 000 5 | 133 x 133 x 52 |
| 96-7220 | Aluminium alloy | Silding Drawer Freezer Racks | 4 x 5 | 564 x 141 x 288.5 | 133 X 133 X 52 |
| 96-7224 | Aluminium alloy | Sliding Drawer Freezer Racks | 4 x 6 | 564 x 141 x 346 | 133 x 133 x 52 |



Frame Type Freezer Racks

- Made of corrosion-resistant stainless steel or aluminium alloy
- Ideal for storage of 1 inch and 2 inch standard size cryogenic boxes
- Fits most standard upright freezers
- Simple structure and lightweight design helps manage cryogenic boxes efficiently
- Can be customized

| Cat # | Material | Description | Capacity (D x H Configuration) | Dimensions (mm) (D x W x H) | Storage box size (mm) |
|---------|-----------------|--------------------------|-----------------------------------|--------------------------------|-----------------------|
| 96-2216 | Aluminium alloy | Frame Type Freezer Racks | 4 × 4 | 552 × 137 × 228 | 133 x 133 x 52 |
| 96-2220 | Aluminium alloy | Frame Type Freezer Racks | 4 x 5 | 552 × 137 × 284 | 133 x 133 x 52 |
| 96-2224 | Aluminium alloy | Frame Type Freezer Racks | 4 x 6 | 552 × 137 × 340 | 133 x 133 x 52 |
| 99-2215 | Stainless steel | Frame Type Freezer Racks | 3 x 5 | 422 × 137 × 284 | 133 x 133 x 52 |
| 99-2216 | Stainless steel | Frame Type Freezer Racks | 4 x 4 | 560 x 137 x 229 | 133 x 133 x 52 |
| 99-2220 | Stainless steel | Frame Type Freezer Racks | 4 x 5 | 560 x 137 x 284 | 133 x 133 x 52 |
| 99-2224 | Stainless steel | Frame Type Freezer Racks | 4 x 6 | 560 x 137 x 339 | 133 x 133 x 52 |
| 99-2225 | Stainless steel | Frame Type Freezer Racks | 5 x 5 | 680 x 137 x 284 | 133 x 133 x 52 |
| 99-2230 | Stainless steel | Frame Type Freezer Racks | 5 x 6 | 696 x 137 x 342 | 133 x 133 x 52 |



Side Access Drawer Freezer Racks

- Made of corrosion-resistant stainless steel
- Ideal for storage of 1 inch and 2 inch standard size cryogenic boxes
- Fits most standard upright freezers
- Can be customized


| Cat # | Material | Description | Capacity (D x H Configuration) | Dimensions(mm) (D x W x H) | Storage box size (mm) |
|---------|-----------------|-------------------------------------|-----------------------------------|-------------------------------|-----------------------|
| 99-5216 | Stainless steel | Side Access Drawer Freezer Racks | 4 x 4 | 560 x 140 x 240 | 133 x 133 x 52 |
| 99-5220 | Stainless steel | Side Access Drawer Freezer Racks | 4 x 5 | 560 x 140 x 300 | 133 x 133 x 52 |
| 99-5225 | Stainless steel | Side Access Drawer Freezer Racks | 5 x 5 | 685 x 140 x 293 | 133 x 133 x 52 |
| 99-5224 | Stainless steel | Side Access Drawer Freezer Racks | 4 x 6 | 560 x 140 x 360 | 133 x 133 x 52 |



Vertical Type Freezer Racks

- Made of corrosion-resistant stainless steel or aluminium alloy
- Ideal for storage of 1 inch and 2 inch standard size cryogenic boxes
- Fits most standard upright freezers
- Especially suitable for liquid nitrogen freezing
- Can be customized

| Cat # | Material | Description | Capacity (D x H Configuration) | Dimensions (mm) (D x W x H) | Storage box size (mm) |
|----------------|-----------------------|---------------|-----------------------------------|--------------------------------|-----------------------|
| 99-62XX | Stainless steel | Vertical Type | 1 x X | * x 140 x 140 | 133 x 133 x 52 |
| 96-62XX | Aluminium alloy | Vertical Type | 1 x X | * x 139 x 140 | 133 x 133 x 52 |
| X (Capacity /L | _) 06, 07, 08, 09, 10 |), 11, 12, 13 | * (Depth /mm) 33 | | 8, 674, 730 |



Sliding Side Drawer Freezer Racks

- Made of high-performance aluminum alloy
- Light weight, convenient operation, and elegant appearance
- Ideal for storage of 1 inch and 2 inch cryogenic boxes and standard freezers
- Handle design enables easy removal from the shelf
- Flanging design reduces scratching
- Can be customized

| Cat # | Material | Description (D | Capacity x H Configuration) | Dimensions (mm) (D x W x H) | Storage box size (mm) |
|---------|-----------------|----------------------------------|--------------------------------|--------------------------------|-----------------------|
| 96-5216 | Aluminium alloy | Sliding Side Drawer Freezer Rack | s 4 x 4 | 564 × 141 × 240.2 | 133 x 133 x 52 |
| 96-5220 | Aluminium alloy | Sliding Side Drawer Freezer Rack | s 4 x 5 | 564 × 141 × 300 | 133 x 133 x 52 |

8-Channel Screw Cap Capper/Decapper & Scanners



8-Channel Screw Cap Capper/Decapper

- Easy-to-use design: Caps or decaps 8 tubes in an individual line Convenient operation: Improves processing throughput of manual procedures with automated column-wise decapping and capping in just three seconds
- Optimal torque level: Ensures reproducibility by capping/ decapping all tubes to the optimal torque level
- Ergonomic design: Allows anyone in the lab to cap or decap tubes with a light weight, easy—to—use and ergonomic design
- Extra-long standby time: Enhances the convenience of actual operation
- Low noise: Enables the capping /decapping noise less than 5 decibels

| Cat# | Lithium battery capacity | Continual operation | Standby time | 8-Channel screw cap capping/Decapping | External dimension(L*W*H) |
|----------|--------------------------|---------------------|--------------|---------------------------------------|---------------------------|
| CKS-1008 | 5500MA | 8H | 30 Davs | 20 | 95mm*46mm*255mm |



CryoKING Handheld Scanner

- Infrared scanning individual cryogenic vials to identify barcodes and 2D barcodes accurately
- Ergonomic handheld design
- Compatible with CryoKING Crown biobanking management software
- Collect sample information rapidly and accurately



CryoKING Scanner-3.0

- Scans one rack in less than 1 second
- Advanced error handling options
- Operates down to -20°C
- Compatible with most cryogenic boxes of major brands
- Project based workflow for easy rack type changes
- Detects missing tubes
- Advanced error handling options

| Cat # | Description | Scan Type | |
|----------|----------------------------|-----------------|--|
| CKS-0101 | Wired Handheld Scanner Gun | Individual Vial | |
| CKS-1301 | AC 100-240V/50Hz/60Hz | Entire Box | |

CryoKING Tube Handling System






>> Features

- Re-array, reformat or cherry pick tubes
- Weigh tubes (tare & gross weights)
- Detect volume with the LevelCheck
- 2D Barcode decoding
- LIMS integration ready
- Work list based operation
- Efficient for both small and large projects
- Label printer and applicator option
- Designed for continuous use and unattended operation
- Encoder based robotics and precision bearings provide the highest reliability
- Integration software toolkits available

Ordering Information

- Electrical: 110-220 VAC 50/60Hz
- Operating Environment: 10° C to 40° C, $10-90^{\circ}$ RH [Standard]; -4° C to 40° C, $10-90^{\circ}$ RH [Cold Room Model]
- System Requirements: Windows 10, 8, 7 2 GB RAM Recommended USB 2.0 ports



Weigh Tubes





Label Printing and Application

| Decode 2D barcodes |
|--|
| 9900 |
| |
| |
| The state of the s |


| Cat # | Capacity | Modules | Dimensions (mm) | Weigjt |
|----------|--------------|--|-------------------|-------------------|
| CKS-1520 | 20 SBS Racks | Standard: Re-array, reformat tubes; Decode 2D barcodes, weigh tubes | 72 x 53 x 40 cm | 27 kg |
| CKS-152X | 20 SBS Racks | Optional: Level check, DMSO hydration check, label, enclosure | Depend on modules | Depend on modules |

Biobanking Management Software

The CryoKING Crown Biobanking Management Software integrates everything from sample management and storage to packaging and transport into one comprehensive management program. It is the tool of choice for streamlined scientific research and biobanking management.

>> Features


- Can be accessed through the CryoKING Cloud service, or hosted locally for use on a laboratory's own LAN
- Designed for use with industry-standard equipment and procedures
- Easy-to-use interface allows users to upload and import existing sample data
- Compatible with CryoKING Scanner-1.0 and Scanner-3.0 for entire box of individual 2D-barcoded cryogenic vials
- Currently available in English and Simplified Chinese (more languages coming soon)






1. Our products (including freezer racks, cryogenic boxes, and cryogenic vials, with more to come) are designed with to be compatible with standard-size equipment and to be suitable for the storage of many types of samples.

2. The Crown visual interface displays a visual representation of every level of a lab's sample storage from freezers to individual samples, allowing users to review their sample storage status at a glance.



Icons and templates can be customized to suit each lab's needs, and custom images can be uploaded for use as cons.



 Transition seamlessly from your previous system b importing your data into the Crown system using Excel c CSV files.

Low Temperature Freezers



Suitable for blood banks, hospitals, sanitation and antiepidemic stations and clinics, electronic industries, university laboratories, military industries, pelagic fishery companies, and drugstores.

-164°C Ecryo Freezer -152°C /-135°C /-105°C Cryo Freezer -86°C Ultra Low Freezer

-86°C Ultra Low Freezer -40°C Ultra Low Freezer -25°C Medical Freezer

Combined Refrigerator &Freezer 2 ~ 8°C Medical Refrigerator 4°C Blood Bank Refrigerator



Liquid Nitrogen Tanks



Made of high quality stainless steel or durable aluminum; intelligent monitoring system and multi alarms to ensure the safety of samples. Different types can be chosen according to different storage requirements.

Filling Series Transport Series Air Transport Series

Special Transport Series Static Storage Series Biobank Series

Portable Dewar Series Laboratory Series



Cell Culture Series

Cell culture involves sampling cells from an organism, placing these cells in an in vitro environment, and observing the growth, reproduction, aging or other processes of the cells. Cell culture has become an important practice in life science research, and plays a central role in the development of the biotechnology industry.

The Biologix Cell Culture Series offers a range of safe and reliable basic supplies for cell transfer, filtering, and other applications.



Cell Culture Flasks







>> Features


- Made of high quality medical grade polystyrene
- Available flask sizes include surface areas of 25 cm², 75 cm², or 175 cm²
- Plug cap and filter cap options available for all flask sizes
- Plug caps completely isolate the culture from the outside environment
- Filter caps prevent cross-contamination while meeting the requirements for gas exchange
- Short, wide, angled neck design for easy access
- Gamma radiation sterilized
- Non-pyrogenic, non-cytotoxic
- DNase & RNase free, human DNA free

Suitable for open culture conditions, effectively preventing cross-contamination while meeting therequirements for gas exchange in cell and tissue culture











Dimensions = mm

| Cat # | Сар Туре | Culture Area | Working Volume | Total Volume | Surface Treatment | Sterile | Unit | Total |
|---------|----------|--------------------|----------------|--------------|-------------------|---------|----------------------|-----------|
| 07-8025 | Filter | 25cm ² | 7ml | 60 ml | Yes | Yes | 5 /Bag, 40 Bags/Case | 200 /Case |
| 07-8075 | Filter | 75cm² | 25ml | 250ml | Yes | Yes | 5 /Bag, 20 Bags/Case | 100 /Case |
| 07-8175 | Filter | 175cm ² | 50ml | 650ml | Yes | Yes | 5 /Bag, 8 Bags/Case | 40 /Case |
| 07-9025 | Plug | 25cm ² | 7ml | 60ml | Yes | Yes | 5 /Bag, 40 Bags/Case | 200 /Case |
| 07-9075 | Plug | 75cm ² | 25ml | 250ml | Yes | Yes | 5 /Bag, 20 Bags/Case | 100 /Case |
| 07-9175 | Plug | 175cm ² | 50ml | 650ml | Yes | Yes | 5 /Bag, 8 Bags/Case | 40 /Case |


Cell Culture Dishes








- Made of high quality medical grade polystyrene
- Lid design promotes effective gas exchange
- External grip for easier handling (07-3035 and 07-3060)
- Raised edges for convenient stacking
- Vacuum plasma surface treatment ensures excellent cell adhesion
- Smooth, transparent materials facilitate observation of cells through a microscope
- Electron beam sterilized
- Non-pyrogenic, non-cytotoxic
- DNase & RNase free, human DNA free







Dimensions = mm


| Cat # | Dish style DxH(mm) | Internal Dimensions D x H (mm) | Culture Area | Working Volume | External Grip | Sterile | Unit |
|---------|--------------------|-----------------------------------|-----------------------|----------------|---------------|---------|--------------------|
| | | | | | | | |
| 07-3035 | 35.00 x 10.00 | 34.30 x 9.30 | 9.40cm ² | 3.00ml | Yes | Yes | 10 /Bag, 500 /Case |
| 07-3060 | 60.00 x 15.00 | 52.80 x 12.80 | 21.50cm ² | 5.00ml | Yes | Yes | 10 /Bag, 500 /Case |
| 07-3100 | 90.00 x 20.00 | 86.26 x 17.70 | 57.50cm ² | 12.50ml | No | Yes | 10 /Bag, 200 /Case |
| 07-3150 | 150.00 x 25.00 | 140.20 x 23.30 | 148.00cm ² | 35.00ml | No | Yes | 10 /Bag, 120 /Case |


Cell Culture Plates





- Constructed of prime virgin polystyrene
- Plate cover design minimizes the risk of cross contamination
- Thin-walled well bottoms increase thermal consistency, reducing the edge effect
- Flat well bottoms provide excellent clarity for microscopic examination of samples
- Raised well rim design around each well enables secure sealing and safeguards against evaporation
- Gamma radiation sterilized
- Non-pyrogenic, non-cytotoxic
- DNase & RNase free, human DNA free






Dimensions = mm

| Cat # | Wells | External Dimensions W x L x H (mm) | Well Dimensions D x H (mm) | Culture Area | Working Volume | Sterile | Unit |
|---------|---------|---------------------------------------|-------------------------------|---------------------|----------------|---------|----------------------|
| 07-6006 | 6-well | 85.40 x 127.60 x 20.20 | 35.00 x 17.50 | 9.6cm² | 3ml | Yes | 1 /Bag, 50 Bags/Case |
| 07-6012 | 12-well | 85.40 x 127.60 x 20.20 | 21.90 x 17.50 | 3.8cm ² | 2ml | Yes | 1 /Bag, 50 Bags/Case |
| 07-6024 | 24-well | 85.40 x 127.60 x 20.20 | 15.50 x 17.50 | 1.9cm ² | 1ml | Yes | 1 /Bag, 50 Bags/Case |
| 07-6048 | 48-well | 85.40 x 127.60 x 20.20 | 9.75 x 17.50 | 0.75cm ² | 0.5ml | Yes | 1 /Bag, 50 Bags/Case |
| 07-6096 | 96-well | 85.40 x 127.60 x 14.40 | 6.50 x 10.80 | 0.33cm ² | 0.2ml | Yes | 1 /Bag, 50 Bags/Case |


Serological Pipettes








- Constructed of highly transparent prime virgin polystyrene
- Available in 1ml, 2ml, 5ml, 10ml, 25ml, 50ml volumes
- Each volume labeled with a different color for ease of use
- Ascending and descending graduations
- Filter at the top of each pipette prevents damage to pipettors
- Electron beam sterilized
- Non-pyrogenic, non-cytotoxic
- DNase & RNase free, human DNA free







 $\mathsf{Dimensions} = \mathsf{mm}$


| Cat # | Working Volume | Color | Sterile | Unit | Total |
|---------|----------------|--------|---------|-------------------------------------|-----------|
| 07-5001 | 1.00ml | Yellow | Yes | 1 /Bag, 200 Bags/Pack, 4 Packs/Case | 800 /Case |
| 07-5002 | 2.00ml | Green | Yes | 1 /Bag, 150 Bags/Pack, 4 Packs/Case | 600 /Case |
| 07-5005 | 5.00ml | Blue | Yes | 1 /Bag, 100 Bags/Pack, 4 Packs/Case | 400 /Case |
| 07-5010 | 10.00ml | Orange | Yes | 1 /Bag, 100 Bags/Pack, 4 Packs/Case | 400 /Case |
| 07-5025 | 25.00ml | Red | Yes | 1 /Bag, 50 Bags/Pack, 4 Packs/Case | 200 /Case |
| 07-5050 | 50.00ml | Purple | Yes | 1 /Bag, 40 Bags/Pack, 4 Packs/Case | 160 /Case |

Cell Culture Slides





- Easy-to-open, highly transparent slide chambers
- Chamber numbering for convenient sample identification
- Polystyrene chamber piece available with 1, 2, 4, or 8 wells
- Packing trays can be used as incubation racks in CO2 incubators
- Slide and chamber apparatus can be easily disassembled without damaging samples
- EO Sterilized




| Cat # | Material Chamber/Slide/Holder | Wells | Color | Growth Area | Working | Sterile | Unit | Total |
|---------|----------------------------------|--------|-------|---------------------|-------------|---------|-----------------------|----------|
| | | | | | | | | |
| 07-2101 | PS /Glass/PP | 1-well | Clear | 9.40cm ² | 2.50-5.50ml | Yes | 6 /Pack, 2 Packs/Case | 12 /Case |
| 07-2102 | PS /Glass/PP | 2-well | Clear | 4.55cm ² | 1.20-2.50ml | Yes | 6 /Pack, 2 Packs/Case | 12 /Case |
| 07-2104 | PS /Glass/PP | 4-well | Clear | 2.13cm ² | 0.50-1.30ml | Yes | 6 /Pack, 2 Packs/Case | 12 /Case |
| 07-2108 | PS /Glass/PP | 8-well | Clear | 0.98cm ² | 0.20-0.60ml | Yes | 6 /Pack, 2 Packs/Case | 12 /Case |


Cell Strainers






- Made of sturdy polypropylene and strong nylon mesh
- Three different pore sizes available: 40µm, 70µm, or 100µm
- Strainers are color-coded according to pore size for convenience
- Fits perfectly into Biologix 50ml Centrifuge Tubes
- Gamma radiation sterilized, non-pyrogenic







 ${\sf Dimensions} = {\sf mm}$

| Cat # | Size | Color | Sterile | Package | Unit | Total |
|---------|-------|--------|---------|------------|-----------------------|-----------|
| 15-1040 | 40µm | Purple | Yes | Individual | 1 /Bag, 100 Bags/Case | 100 /Case |
| 15-1070 | 70μm | Orange | Yes | Individual | 1 /Bag, 100 Bags/Case | 100 /Case |
| 15-1100 | 100µm | Green | Yes | Individual | 1 /Bag, 100 Bags/Case | 100 /Case |

Sterile Syringe Filters & Vacuum Units





- \bullet Syringe filters are suitable for 1ml~150ml liquid or solution, vacuum units are suitable for 150ml~2L liquid or solution
- Single Use and Gamma Ray pre-sterilized
- Sterile graded 0.2um membrane: LRV > 7 for Brevundimonas diminuta, ATCC19146.
- Non-pyrogenic, low extractables and surfactant-free ensures purity of filtrate.
- Extreme-low Protein Binding <1.5% (PES)
- Excellent Throughput, for faster filtration and more process volume.

| Cat # | Description | Unit |
|-------------|--|-------------------------|
| 08-PES1322R | Syringe filters. Membrane: PES, pore size: 0.22µm, diameter:13mm. Female luer lock + male luer slip. Extreme-low protein binding. Individually packed | 100 /Pack, 8 Packs/Case |
| 08-PES2522R | Syringe filters. Membrane: PES, pore size: 0.22µm, diameter:25mm. Female luer lock + male luer slip. Extreme-low protein binding. Individually packed | 100 /Pack, 4 Packs/Case |
| 08-PES3022R | Syringe filters. Membrane: PES, pore size: 0.22µm, diameter:30mm. Female luer lock + male luer slip. Extreme-low protein binding. Individually packed | 100 /Pack, 4 Packs/Case |
| 08-PES3322S | Syringe filters. Membrane: PES, pore size: 0.22µm, diameter:33mm. Double luer lock. Extreme-low protein binding. Individually packed | 100 /Pack, 4 Packs/Case |
| 08-PES3345S | Syringe filters. Membrane: PES, pore size: 0.45µm, diameter:33mm. Double luer lock. Extreme-low protein binding. Individually packed | 100 /Pack, 4 Packs/Case |
| 08-PES3322M | Syringe filters. Membrane: PES, pore size: 0.22µm, diameter:33mm. Female luer lock + male luer slip. High pressure resistance. Individually packed | 100 /Pack, 4 Packs/Case |
| 08-PES3345M | Syringe filters. Membrane: PES, pore size: 0.45µm, diameter:33mm. Female luer lock + male luer slip. High pressure resistance. Individually packed | 100 /Pack, 4 Packs/Case |
| 08-CA3022R | Syringe filters. Membrane: CA, pore size: 0.22µm, diameter:30mm. Female luer lock + male luer slip. High biocompatibility. Individually packed | 100 /Pack, 4 Packs/Case |
| 08-CA3045R | Syringe filters. Membrane: CA, pore size: 0.45µm, diameter:30mm. Female luer lock + male luer slip. High biocompatibility. Individually packed | 100 /Pack, 4 Packs/Case |
| 09-PES2250 | Vacuum filtration unit. Membrane: PES, pore size: 0.22µm, diameter:90mm, funnel volume: 500ml, receiver volume: 500ml | 12 /Case |
| 09-PES2250F | Vacuum filtration unit-funnel. Membrane: PES, pore size: $0.22\mu m$, diameter: 90mm, funnel volume: 500ml | 24 /Case |
| 09-CA2250 | Vacuum filtration unit. Membrane: CA, pore size: 0.22µm, diameter:90mm, funnel volume: 500ml, receiver volume: 500ml | 12 /Case |
| 09-CA2250F | Vacuum filtration unit-funnel. Membrane: CA, pore size: 0.22μm, diameter:90mm, funnel volume: 500ml. | 24 /Case |

Cell Scrapers & Lifters






>> Features


Cell Scrapers

- Made of polyethylene
- Available in two sizes: 18cm or 25cm in length,
- Pivoting blade makes precise cutting and scraping movements
- Gamma radiation sterilized, non-pyrogenic

Cell Lifters

- Made of polyethylene
- One-piece design and chiseled blade provide improved control while minimizing damage to cells
- Gamma radiation sterilized, non-pyrogenic





| Cat # | Name | Handle Length | Sterile | Package | Unit | Total |
|---------|--------------|---------------|---------|------------|-----------------------|-----------|
| 70-2180 | Cell Lifter | 18cm | Yes | Individual | 1 /Bag, 100 Bags/Case | 100 /Case |
| 70-1180 | Cell Scraper | 18cm | Yes | Individual | 1 /Bag, 100 Bags/Case | 100 /Case |
| 70-1250 | Cell Scraper | 24.3cm | Yes | Individual | 1 /Bag, 200 Bags/Case | 200 /Case |

Dimensions = mm



PCR Tubes






>> Features

- Made of high quality medical grade polypropylene
- Ultra thin-wall tubes provide faster heating and cooling
- High degree of transparency allows for easy observation of samples
- Available with domed caps or flat caps
- Domed caps ensure efficient heat transfer
- Flat caps are ideal for Real-Time PCR (qPCR)
- Unique cap design prevents opening during cycling and incubation stages
- DNase, RNase, and endotoxin free

- Fits most standard thermocyclers
- Can withstand centrifugation of up to 6,000 RCF




Dimensions = mm

| Cat # | Volume | Color | Сар | Sterile | Unit | Total |
|----------|--------|--------|-------------|---------|--------------------------------|-------------|
| 60-0082 | 0.2ml | Clear | Flat-top | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0182* | 0.2ml | Pink | Flat-top | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0282* | 0.2ml | Green | Flat-top | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0382* | 0.2ml | Blue | Flat-top | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0482* | 0.2ml | Yellow | Flat-top | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0002 | 0.2ml | Clear | Dome-shaped | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0102* | 0.2ml | Pink | Dome-shaped | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0202* | 0.2ml | Green | Dome-shaped | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0302* | 0.2ml | Blue | Dome-shaped | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |
| 60-0402* | 0.2m | Yellow | Dome-shaped | No | 1,000 /Bag/Pack, 10 Packs/Case | 10000 /Case |

^{*} Call for details.

8-Strip PCR Tubes






>> Features

- Made of high quality medical grade polypropylene
- Ultra thin-wall tubes provide faster heating and cooling
- High degree of transparency allows for easy observation of samples
- Available with domed caps or flat caps
- Domed caps ensure efficient heat transfer
- Flat caps are ideal for Real-Time PCR (qPCR)
- Unique cap design prevents opening during cycling and incubation stages
- DNase, RNase, and endotoxin free

- Fits most standard thermocyclers
- Can withstand centrifugation of up to 6,000 RCF
- Packaged in resealable bags




 ${\sf Dimensions} = {\sf mm}$

| Cat# | Volume | Color | Сар Туре | Sterile | Unit |
|----------|--------|--------|----------|---------|--|
| 60-0088 | 0.2ml | Clear | Flat | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0188* | 0.2ml | Pink | Flat | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0288* | 0.2ml | Green | Flat | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0388* | 0.2ml | Blue | Flat | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0488* | 0.2ml | Yellow | Flat | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0008 | 0.2ml | Clear | Domed | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0108* | 0.2ml | Pink | Domed | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0208* | 0.2ml | Green | Domed | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0308* | 0.2ml | Blue | Domed | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |
| 60-0408* | 0.2ml | Yellow | Domed | No | 125 Strip Tubes+125 Strip Caps/Bag/Pack, 10 Packs/Case; 1250 sets/Case |

^{*} Call for details.



PCR Plates






>> Features

- Made of high quality medical grade polypropylene
- Ultra-thin tube walls with uniform thickness for accelerated heat transfer, shorter cycle times, and more consistent results
- Alphanumeric well labeling for convenient identification
- DNase, RNase, and endotoxin free

- Fits most standard thermocyclers
- Can withstand centrifugation of up to 6,000 RCF



Half-skirted plates are easy to insert and remove



Dimensions = mm

| Cat # | Well | Skirt | Color | Sterile | Unit | Total |
|----------|---------|--------------|--------|---------|------------------------------------|-----------|
| 60-0056 | 96-well | Half-Skirted | Clear | No | 5 /Bag, 5 Bags /Pack, 4 Packs/Case | 100 /Case |
| 60-0156* | 96-well | Half-Skirted | Pink | No | 5 /Bag, 5 Bags /Pack, 4 Packs/Case | 100 /Case |
| 60-0256* | 96-well | Half-Skirted | Green | No | 5 /Bag, 5 Bags /Pack, 4 Packs/Case | 100 /Case |
| 60-0356* | 96-well | Half-Skirted | Blue | No | 5 /Bag, 5 Bags /Pack, 4 Packs/Case | 100 /Case |
| 60-0456* | 96-well | Half-Skirted | Yellow | No | 5 /Bag, 5 Bags /Pack, 4 Packs/Case | 100 /Case |

^{*} Call for details.

PCR Plate Compatibility

| Biolo | ogix Plates 🕨 | 60-0056 (Half-Skirted) |
|-----------------------|------------------|---------------------------|
| Brand | Model | |
| Real Time PCR Thermal | Cyclers 7000 | 1/ |
| Applied Biosystems | 7300 | |
| | 7500 | 1/ |
| | 7700 | 1/ |
| | 7900HT | ${}$ |
| | StepOne | V |
| BioRad | iCycler | |
| 2.0.1.0.0 | MyiQ | • |
| | iQ5 | |
| Eppendorf MasterCyc | cler ep Realplex | • |
| MJ Research | Opticon | |
| | Opticon 2 | |
| | Chromo 4 | |
| | MX-4000 | |
| | MX-3000 | |
| Techne | Quantica | |
| PCR Thermal Cycl | | |
| Applied Biosystems | Veriti | |
| | 9600 | |
| 9700 | | |
| 7800FAST | | |
| 2720 | | |
| Biometra | Uno | |
| Uno II | | |
| T1 Thermal Cycler | | |
| Tgrasient | | |
| Trobot | | |
| BioRad | iCycler | |
| MyCycler | | |
| Corbett Research | PalmCycler 96 | |
| PalmCycler 384 | | |
| Eppendorf | Mastercycler | |
| Mastercycler ep | | |
| Mastercycler M384 | | |
| Flexi | Gene | V |
| Genius | | |

| | Biologix Plates > | 60-0056 (Half-Skirted) |
|------------------|-----------------------|---------------------------|
| Brand | Model | |
| Ericom | Single Block | |
| | Twin Block | |
| | Deltacycler I | |
| ThermoHybaid | PCR Sprint | |
| | PCR Express | |
| | Multiblock System | |
| | Touchdown | |
| | Omnigene | |
| | Omn-E | |
| | рхе | |
| | рх2 | |
| MJ Research | PTC-100 | |
| | PTC-200 | |
| | PTC-225 Tetrad | $\sqrt{}$ |
| | Dyad/Disciple | $\sqrt{}$ |
| MWG | Primus 96 | $\sqrt{}$ |
| | Primus 384 | |
| Stratagene | Robocycler | |
| TaKaRa | TP240 | |
| | TP3000 | |
| Techne Toucho | gene Gradient (TC512) | $\sqrt{}$ |
| | Genius (TC412) | V |
| | Touchgene X | |
| Sequencers | | |
| Applied Biosyste | ms 3100 | |
| 3130 | | |
| 3700 | | |
| 3730/3730x | | |

This chart is intended only as a general guide. Compatibility may vary based on your particular configuration and features. For further help with product selection, please contact us.

Microbiology Series




Petri Dishes






>> Features

- Made of medical grade polystyrene
- Highly transparent for microscopic analysis
- Three sizes available: 60mm, 90mm, or 150mm in diameter
- 66-1701 fit automatic systems
- Vented lid for better gas exchange
- Dishes are stackable
- EO Sterile






Dimensions = mm




Inoculating Loops & Needles







- Made of medical grade polystyrene
- Loops available in 1µL and 10µL sizes
- The opposite end of each loop can be used as a needle
- Needles are suitable for removal of single-colony specimens
- Gamma radiation sterilized, non-pyrogenic
- Packaged in resealable bags




| Cat # | Size | Description | Color | Sterile | Unit | Total |
|---------|------|-------------|--------|---------|-----------------------|------------|
| 65-0001 | 1µl | Loops | Blue | Yes | 25 /Bag, 40 Bags/Case | 1000 /Case |
| 65-0010 | 10µl | Loops | Orange | Yes | 25 /Bag, 40 Bags/Case | 1000 /Case |
| 65-0002 | - | Needles | Yellow | Yes | 25 /Bag, 40 Bags/Case | 1000 /Case |

Cell Spreaders





- Polystyrene
- Smooth surface quality and angled end reduce the risk of damaging the culture medium
- L-shape allows application on the entire Culture Dish or plate surface
- Two package options: individually wrapped (sterile), or 10 per bag (sterile, resealable)
- Gamma radiation sterilized



 ${\sf Dimensions} = {\sf mm}$

| Cat # | Color | Sterile | Unit | Total |
|---------|--------|---------|-----------------------|-----------|
| 65-1001 | Orange | Yes | 1 /Bag, 500 Bags/Case | 500 /Case |
| 65-1010 | Orange | Yes | 10 /Bag, 50 Bags/Case | 500 /Case |




Drosophila Vials





- Made of premium grade polystyrene
- Glass-like clarity for sample observation
- Disposable vials are a safe and cost-effective alternative to glass vials
- Narrow or wide vials available
- Both narrow and wide vials are available in bulk or tray packaging
- Tray packs are great for quick, effortless filling of vials



Dimensions = mm

| Cat # | Description | Dimensions (mm) | Sterile | Packaging | Unit | Total |
|---------|-------------|-----------------|---------|-----------|-------------------------|-----------|
| 51-0500 | Narrow | 25 x 95 | No | Bulk | 25 /Bag, 20 Bags/Case | 500 /Case |
| 51-0501 | Narrow | 25 x 95 | No | Tray | 100 /Tray, 5 Trays/Case | 500 /Case |
| 51-0800 | Wide | 28.5 x 95 | No | Bulk | 25 /Bag, 20 Bags/Case | 500 /Case |
| 51-0801 | Wide | 28.5 x 95 | No | Tray | 100 /Tray, 5 Trays/Case | 500 /Case |

Drosophila Bottles






>> Features

- Made of clear polypropylene
- Autoclavable
- Transparent bottle with square bottom
- Firm, stiff neck for easier handling






Dimensions = mm

| Cat # | Sterile | Packaging | Unit | Total |
|----------|---------|-----------|------------------------|-----------|
| 51-17720 | No | Bulk | 200 /Bag, 1 Bag/Case | 200 /Case |
| 51-17721 | No | Tray | 25 /Tray, 8 Trays/Case | 200 /Case |

Drosophila Plugs





- Made of high density foam
- Causes minimal shedding and reduces evaporation
- Maintains tight fit after repeated use and removal
- Fits Biologix Drosophila vials and bottles

Flat surface accepts can be marked













51-17725

| Cat # | Applied # | Fits | Sterile | Unit | Total |
|----------|----------------|-------------------------|---------|------------------------|------------|
| 51-0505 | 51-0500/0501 | Narrow Drosophila Vials | No | 100 /Bag, 10 Bags/Case | 1000 /Case |
| 51-0805 | 51-0800/0801 | Wide Drosophila Vials | No | 100 /Bag, 10 Bags/Case | 1000 /Case |
| 51-17725 | 51-17720/17721 | Drosophila Bottles | No | 50 /Bag, 20 Bags/Case | 1000 /Case |



Sample Collection and Detection Products

Biologix supplies various sample collection and detection tools of complete specifications. Suitable for the collection, transportation, preservation, and detection of human samples and environment samples.





Sample Collection and Detection Products

Disposable Virus Collection Tube & Transportation, Preservation Medium



6ml or 12ml tubes with 2ml or 3ml medium.

All available: Oral or nasal swabs, inactivated or non-inactivated

Package: Racked tubes or individually wrapped kits available.



Sample Collection Tubes



6ml or 10ml tubes without medium, self-standing bottom. Sterile.

Disposable Swabs



Oral or nasal swabs, individually wrapped; Sterile. Breaking point: 30mm or 80mm.

Disposable Saliva Sample Collectors



10ml tubes with inactivated or non-inactivated medium. Biohazard disposable bags available. Individually wrapped.



Liquid Handling



dPette Digital Pipette

Features

- Comprehensive range of liquid handling protocols
- Three aspiration/dispensing speeds
- Motor-drive with built-in error control
- Efficient lithium-ion rechargeable battery
- Autoclavable lower barrel





Micro PettePlus Pipettes

Features

- Convenient one-handed operation
- Fully autoclavable
- Compatible with most standard pipette tips

Pipette Stand (01-2306)

• Holds up to six single- or multi-channel pipettes

| Cat # | Туре | Channels | Volume Range (μL) | Tip Volume |
|---------|--------------------------|----------|-------------------|------------|
| 01-241X | dPette Digital Pipette | 1 | 0.5-10μL | 0.01µL |
| 01-242X | dPette Digital Pipette | 1 | 5-50µL | 0.1µL |
| 01-243X | dPette Digital Pipette | 1 | 30-300μL | 1μL |
| 01-244X | dPette Digital Pipette | 1 | 100-1000μL | 5μL |
| 01-2101 | MicroPette Plus Pipettes | 1 | 0.5-10μL | 0.1µL |
| 01-2104 | MicroPette Plus Pipettes | 1 | 10-100μL | 1µL |
| 01-2102 | MicroPette Plus Pipettes | 1 | 20-200μL | 1µL |
| 01-2110 | MicroPette Plus Pipettes | 1 | 100-1000μL | 5µL |
| 01-2115 | MicroPette Plus Pipettes | 1 | 1000-5000μL | 50μL |
| 01-2112 | MicroPette Plus Pipettes | 1 | 0.1-2.5µL | 0.05µL |
| 01-2122 | MicroPette Plus Pipettes | 1 | 2-20µL | 0.5µL |
| 01-2155 | MicroPette Plus Pipettes | 1 | 5-50µL | 0.5µL |
| 01-2801 | MicroPette Plus Pipettes | 8 | 0.5-10μL | 0.1µL |
| 01-2805 | MicroPette Plus Pipettes | 8 | 5-50µL | 0.5µL |
| 01-2830 | MicroPette Plus Pipettes | 8 | 50-300µL | 5µL |

Liquid Handling



Levo ME Pipette Filler

Features

- Ergonomic shape and intuitive operation
- Pressure-sensitive buttons: increasing pressure increases aspiration/dispensing speed
- Autoclavable hydrophobic filter prevents over-aspiration
- LED light indicates aspiration/dispensing, charging, and low battery

| Cat # | Туре | Aspiration Speeds | Dispensing Speeds | Battery Life | Compatible Pipette Types |
|---------|---------|--|--------------------|------------------------------|--|
| 01-221X | Levo ME | aspirates 25mL of water in < 7 seconds | gravity dispensing | 8+ hours of intermittent use | glass or plastic pipettes (0.1-100mL), Pasteur pipettes |



Levo Plus Pipette Filler

- Six aspiration and dispensing speeds
- LCD display shows speed settings and remaining battery life
- Efficient lithium-ion battery can last for up to 8 hours of intermittent use
- Integral, autoclavable silicone adapter suits use with 0.1-100 mL pipettes

| Cat # | Type | Aspirate Speeds | Dispense Speeds | Battery Life | Compatible Pipette Types |
|---------|-----------|------------------|---------------------------|------------------------------|--|
| 01-220X | Levo Plus | 25ml<5s(8 shift) | Motor (6 shift) / Gravity | 8+ hours of intermittent use | glass or plastic pipettes (0.1-100mL), Pasteur pipettes |

Collection System



SAFEVAC Vacuum Aspiration System

Features

- Integrated pump with adjustable vacuum control
- Self-closing connectors for quick installation and dismantling, prevention of drips, and containment of aerosols
- Liquid level sensor prevents overflow
- Single-channel and multi-channel nozzles and other accessory options available

| Cat # | Vacuum range | Pump flow rate | Dimensions [W x D x H] | Weight | Input voltage | Bottle volume |
|---------|--------------|----------------|---------------------------|--------|-----------------------|---------------|
| 01-270X | 0 - 600 mbar | 1-15 ml/s | 180 x240 x 450mm | 3 kg | 100-240 VAC, 50/60 Hz | 4L |



Laboratory Fluid Collection System

- Filmed inner wall
- Connection type: Can be connected with fluid instruments of various brands
- Funnel type: Equipped with the valve funnel
- High sealing performance
- Equipped with super efficient filer with connector
- Optional: Liquid level alarm



| Cat # | Product | Description | Unit |
|---------|-----------------------------------|---|------------|
| 05-2720 | Connection type | One 20L security bucket with inner film , One inner cover, One outer cover, One set of liquid-phase connector, One set of efficient filter connector, One inorganic or organic filter | 1 Set/Case |
| 05-2721 | Funnel type | One 20L security bucket with inner film , One inner cover, One outer cover, One funnel with a cover, One set of efficient filter connector, One inorganic or organic filter | 1 Set/Case |
| 05-2791 | Super efficient filter | norganic chemical filter | 1 Set/Case |
| 05-2792 | Super efficient filter | Organic chemical filter | 1 Set/Case |
| 05-2793 | Liquid lever alarm | Tactile liquid lever alarm | 1 Set/Case |
| 05-2797 | Liquid carrying cart(standard) | Adjustable size; one cart and one tray | 1 Set/Case |
| 05-2798 | Liquid carrying cart(intelligent) | djustable size; one cart and one tray; equipped with intelligent module | 1 Set/Case |



Classic Magnetic Stirrer

Features

- Powerful rare-earth permanent magnet
- PET stirring plate provides excellent chemical resistance
- Stand included



FlatSpin Magnetic Stirrer

Features

- Ultra-flat, space-saving design
- Automatic rotation-reversing stirring for better mixing results
- Slip-resistant and stable



Magnetic Stirrer Hotplate

Features

- Digital control system
- ABS casing is fire-retardant and resistant to weak acids and bases
- Ceramic-coated stainless steel hotplate is chemically resistant and scratch-proof
- Safety feature: HOT caution indicator is displayed at 50°C and above



EcoStir Square Plate Magnetic Stirrer

- Compact design with small footprint, 12V voltage available for mobile labs.
- Square work plate supports 1.5L stirring capacity
- Holder with clamp for pH electrode or another sensor
- Default PET top provides excellent chemical resistance
- Extra colorful stickers satisfy DIY design by yourself
- Maintenance free brushless DC moto

| Cat # | Туре | Stirring plate size [mm] | Motor input rating [W] | Motor output rating [W] | Max. stirring quantity [L, H ₂ O] | Speed range [rpm] | Overheating Protection [°C] |
|----------|----------|-----------------------------|---------------------------|-------------------------|---|----------------------|--------------------------------|
| 01-330X | Classic | Ø 120 | 7.2 | 1.2 | 1.5 | 300-2000 | - |
| 01-320X | FlatSpin | Ø 90 | 5 | 3 | 0.8 | 15-1500 | - |
| 01-310X | Hotplate | Ø 135 | 5 | 3 | 3 | 100-1500 | 320 |
| 01-330XA | EcoStir | Ø 120 | 7.2 | 1.2 | 1.5 | 300-2000 | 5-40 |

Mixers & Rotators



Vortex Mixer

Features

- Orbital mixing movement
- Touch and continuous operation settings
- Cast aluminum base and specially designed vacuum suction feet keep the mixer stable
- Various adapters and accessories available for use with tubes of different sizes (Φ6mm, Φ10mm, Φ12mm, Φ16mm, Φ20mm test tubes and vessels <99mm in diameter)



Classic Rotators (axial) 01-130X



Classic Rotators (radial) 01-120X

Classic Rotators

- End-over-end rotation
- Compact, space saving design
- Various adapters available for use with 1.5mL-50mL tubes

| Cat # | Туре | Speed range [rpm] | Operation mode | Voltage [V] | Power [W] |
|---------|-----------------|-------------------|----------------------------|-----------------|-----------|
| 01-110X | Vortex Mixer | 0-2500 | Touch operation/Continuous | 110-120/220-230 | 60 |
| 01-120X | Classic Rotator | 0-80 | Continuous | 100-240 | 20 |
| 01-130X | Classic Rotator | 0-80 | Continuous | 100-240 | 20 |

LED Digital Orbital and Linear Shaker



- Electronic adjustment of speed and timer
- Double LEDs for independent display of time and speed
- Speed range of 40-200rpm
- Fixed board face, max. load capacity is 2kg with platform
- Continuous operation (endless)
- Electronic time switching clock: 1min-19h59min
- Wide range of platforms for use with a variety of vessels
- Work plate dimension (W×D): 268x268mm
- Max. load capacity (with platform): 2kg
- Timer display: LED
- Time settings range: 1min-l 9h59min
- Power: 20W
- Motor input power: 20WMotor output power: 16W
- Dimension (W×D×H): 294x348x123mm
- Weight: 3.1kg
- Permissible ambient temperature and humidity: 5-40℃,80%RH

| 01-621XA | 20mm | 40-200rpm | LED | DC motor | IP21 | 100-240V,50/60HZ |
|----------|---------------------|-------------|---------------|------------|------------------|------------------|
| Cat # | Orbital diameter | Speed range | Speed display | Motor type | Protection class | Voltage |

Thermocyclers



Gradient Thermocycler

Features

- Thermal gradient: program different temperatures across different areas of the heating block
- Consistent and uniform sealing across entire heating block
- Intuitive touch-screen user interface



Thermocycler

- Consistent and uniform sealing across entire heating block
- Intuitive touch-screen user interface

| | Туре | Heating Temperature Range [°C] | | | Gradient Temperature Setting Range [°C] | Gradient Range [°C] |
|---------|--------------------------|-----------------------------------|-----|---|--|------------------------|
| 01-700X | Gradient Thermocycler | 4-99 | 0.1 | 3 | 30-99 | 1-42 |
| 01-710X | Thermocycler | 4-99 | 0.1 | 3 | - | - |



Palm Micro Centrifuge

Features

- Quick acceleration and deceleration
- Rotors for microcentrifuge tubes, PCR tubes, and PCR 8-strip tubes included.
- Quiet, low-noise design

| Cat # | Туре | Speed Range [rpm] | Max. RCF [×g] | Rotor Capacity |
|---------|--------------------------|-------------------|---------------|--|
| 01-800X | Palm Micro Centrifuge | 7000 | 2680 | 0.2ml/0.5ml/1.5ml/2ml x 8; 0.2ml PCR tubes x 16 or 0.2ml PCR 8-strips x 2 |



Mini and Mini-Plus High Speed Centrifuges

Features

- Brushless DC motor quickly accelerates rotor to the set speed
- Rotors are designed for 1.5ml and 2ml microcentrifuge tubes; adapters for 0.2ml and 0.5ml tubes are also available

Specifications

01-821X (Mini-Plus) with 01-8291 plastic rotor kit*

• Low noise: 54dB; delivers fast acceleration and deceleration

| Cat # | Туре | Speed Range [rpm] | Max. RCF [×g] | Rotor Capacity | Speed Accuracy [rpm] | Temperature Range [°C] |
|---------|----------------------------|------------------------------|--------------------------|----------------------------------|-------------------------|---------------------------|
| 01-821X | Mini High Speed Centrifuge | 500-15000; increment: 100 | 15100; increment: 100 | 0.2ml/0.5ml/ 1.5ml/2.0ml x 12 | ±20 | N/A |

Centrifuges



High Speed Micro-Centrifuge

Features

- Bio-safe rotor prevents sample overflow and ensures samples do not leak during the centrifuging process
- High strength aluminum alloy rotor is autoclavable and highly durable
- Four different rotors are available, assembled with 1.5ml/2ml x 24

| Cat # | Туре | Speed Range [rpm] | Max. RCF [×g] | Rotor Capacity | Speed Accuracy [rpm] | Temperature Range [°C] |
|---------|--------------------------------|----------------------------|-------------------------|--|-------------------------|---------------------------|
| 01-830X | High Speed Micro-Centrifuge | 200-15000; increment:10 | 21380; increment: 10 | 5ml×18; 0.2ml/0.5ml/ 1.5ml/2ml × 24; 0.5ml × 36;PCR8 × 4 | ±20 | N/A |

High-speed Refrigerated Micro Centrifuge

Features

 Centrifuge chamber door with automatic adsorption door lock design,can be closed with a single finger.

 Precooling, the compressor will automatically start cooling as long as the set temperature is lower than the ambient temperature.

 Sound indicator to prompt that the upper cover is closed to maintain the temperature inside the chamber, after centrifugation.

 Advanced R290 flurine free refrigerant ,low caron and environmental protection.

 Up to 9 program storage capacity according to experimental needs.

- Biosafety rotor, three-layer seal for biosafety.
- High-strength aluminum alloy, three layered biosafety to withstand numerous high temperature and high pressure sterilizations.
- Five rotors with 0.2ml,0.5ml,1.5ml,2ml and 5ml centrifuge tubes/PCR tubes capacity.
- An external condensation tank prevents condensate within the centrifuge chamber.
- Other functions: Rotation speed/RCF switchover, inching operation,running process indicator, sound alert.
 Acceleration/deceleration in 9 steps, program storage capacity.



| Cat # | Туре | Max. Speed | Max. RCF | Rotor | Temperature setting |
|---------|--|---|-------------------------|--|---------------------------------------|
| 01-832X | High-speed Refrigerated Micro Centrifuge | 15,000rpm (200-15,000rpm) in 100rpm increments. | 21,380xg, step: 10xg | 1.5/2mL×24,0.5mL×36, PCR8×4,5mL×12, 5mLx18 | 30s-99 min; HOLD (continuous running) |



Clinical Centrifuge

Features

- Max. speed 4000rpm.
- Rotor speed can be set in RPM/RCF, with 3 programs storage capacity for routine applications.
- Max. capacity 50ml×4, Adapt to different sizes of blood collection tubes and centrifuge tubes suitable for swing out rotor and angle rotor. Automatic internal diagnosis and self-detection.
- Maintenance-free brushless DC motor.
- Noise:≤62dB(A)
- Power:Single phase,100V-240V, 50Hz/60Hz.



Clinical Centrifuge

- Speed range of 300-4000rpm.
- 3 programs storage capacity for routine applications.
- Automatic internal diagnosis.
 Maintenance free brushless DC motor.
- Noise:≤60dB(A)
- Power:Single phase,100V-240V, 50Hz/60Hz

| Cat # | Туре | Speed range | Max. RCF | Rotor | Speed accuracy |
|---------|---------------------|----------------------------------|----------|---|----------------|
| 01-860X | Clinical Centrifuge | 500-4000rpm increment: 100rpm | 2500×g | Swing6×15ml Swing4×50ml Angle24×15ml, Angle12×10ml/8×15ml | ±100rpm |
| 01-870X | Clinical Centrifuge | 300-4000rpm increment: 100rpm | 1900×g | Angel8x15ml(Default), Angel10mlx12 | ±100rpm |

Dry Bath Incubators



Dry Bath Incubator

Features

- A broad temperature range, from room temperature to 120°C
- Continuous and timed operation settings
- Over-temperature protection
- Various interchangeable heating blocks are available (0.2ml x 54, 0.5ml x 40, 1.5/2ml x 40, 5/15ml x 28, 50ml x 8, 0.2ml+0.5ml+1.5/2ml, 96 or 384 well plate)

| Cat # | Туре | Functions | Temperature uniformity [°C] | Overheating protection | Max heating rate [°C/min] |
|---------|-----------------------|-----------|-----------------------------|------------------------|------------------------------|
| 01-910X | Dry Bath Incubator | Heating | Max ± 0.5 | 140 ℃ | 5.5 |



Hot/Cold Dry Bath Incubator

- Excellent temperature control system
- Can maintain hot or cold temperatures
- Supports creating and saving custom operation programs
- Various interchangeable heating blocks are available (24 x 0.5ml, 24 x 1.5ml, 24 x 2ml, 8 x 5ml, 8 x 15ml, 4 x 50ml, well plate, 0.2 x 96 PCR plate)

| Cat # | Туре | Functions | Temperature range | Temperature control accuracy [°C, at 20–45°C] | Max heating rate [°C/min] |
|---------|----------------------|---------------------|---|---|------------------------------|
| 01-920X | Hot/Cold Dry Bath | Heating & cooling (| Heating:room temp.~110°C Cooling:room temp.~ below room temp. 25°C | ± 0.5 | 5.5 |

Dry Bath Incubators



Shaker Incubator

Features

- Excellent temperature control system
- Can maintain hot or cold temperatures
- Orbital shaking movement function
- Supports creating and saving custom programs
- Various interchangeable heating blocks are available (24 x 0.5ml, 24 x 1.5ml, 24 x 2ml, 8 x 5ml, 8 x 15ml, 4 x 50ml, well plate, 0.2 x 96 PCR plate)

| Cat # | Туре | Functions | Temperature range | Max heating rate [°C/min] | Max cooling rate [°C/min] | Shaking speed range [rpm] |
|---------|---------------------|-------------------------------|--|------------------------------|--|---------------------------------|
| 01-930X | Shaker Incubator | Heating, cooling, & mixing | Heating:room temp.~100°C Cooling:room temp.~ below room temp. 25°C | 5.5 | 5 (110°C - room temp.) 0.5 (below room temp.) | 300–1500 |





Elite Series Dry Bath Incubators

- PID temperature control maintains temperature accuracy
- Leak-proof molded aluminum alloy chamber with PTFE coating
- Multi-functional: can also be used as a water bath or bead bath
- Various interchangeable heating blocks are available (0.2ml tubes x 64 or PCR 8-strip tubes x 8, 0.5ml tubes x 20, 1.5ml or 2.0ml tubes x 20, 15ml x centrifuge tubes x12, 50ml x centrifuge tubes x4)
- * Blocks are not included, need to buy separatedly.

| Cat # | Туре | Display | Temperature Control Range [℃] | Temperature Increment [℃] | TemperatureAccuracy at 37 [℃] | Operation Modes |
|---------|---|-------------|----------------------------------|------------------------------|-------------------------------|---------------------|
| 03-4221 | Elite Series One-Block Dry Bath Incubator | LCD Display | 5 above ambient to | 0.10 | . 20 | Continuous or timed |
| 03-4222 | Elite Series Dual Block Dry Bath Incubator | LCD Display | 150 | 0.10 | ±20 | up to 99 hr:59 min |

Dry Bath Incubators

Dry bath

Features

- Overheating protection.
- USB interface can be used to trace historical data.
- Sound reminder function.
- External temperature sensor.
- Block equipped with a lid for heat preservation and prevent pollution.
- Knob adjustment is easy to operate.



| Cat # | Туре | Temperature setting range[°C] | Output Current / Inc [mA] | Power [w] | Operating temperature[°C] |
|---------|----------|----------------------------------|---------------------------|-----------|---------------------------|
| 01-900X | Dry bath | 25~105 | 25-90:±0.3 90-105:±0.6 | 100 | +10~40 |
| 01-9A0X | Dry bath | 25~105 | 25-90:±0.3 90-105:±0.6 | 200 | +10~40 |
| 01-950X | Dry bath | 25~150 | 25-90:±0.3 90-105:±0.6 | 100 | +10~40 |
| 01-960X | Dry bath | 25~150 | 25-90:±0.3 90-105:±0.6 | 200 | +10~40 |




Mini Dry bath

- Light in weight
- LCD display of both temperature and time.
- Rapid calibration support
- Overheating protection
- Optional blocks with different capacities are available
- Safe and stable
- USB interface to save the data
- Mini 01-970x with hot lid provided to preserve the generated heat

| Cat # | Туре | Temperature setting range[°C] | Temperature control accuracy[°C] | Power [w] | Operating temperature[°C] |
|---------|---------------|----------------------------------|-------------------------------------|-----------|---------------------------|
| 01-970X | Mini Dry bath | 0~100 | ±0.5 | 60 | +10~40 |
| 01-980X | Mini Dry bath | 0~100 | ±0.5 | 60 | +10~40 |
| 01-990X | Mini Dry bath | 25~150 | ±0.5 | 60 | +10~40 |

Electrophoresis








300V Power Supply


Power Supplies

Features

- Stackable and programmable
- Multiple safety features: no-load detection, shrouded plugs and sockets, leakage detection, over-temperature protection, crossover protection, full floating, over-current/over-voltage detection

| Cat # | Туре | Output Voltage / Inc [V] | Output Current / Inc [mA] | Max Power [W] |
|---------|------------------------|--------------------------|---------------------------|---------------|
| 03-1300 | Mini 300V Power Supply | 10-300/1 | 10-400/1 | 60 |
| 03-2310 | 300V Power Supply | 5-300/1 | 1-700/1 | 150 |





Electrophoresis Tanks

- Injection molded, leak-proof tank
- No tapes, clamps, or springs required
- Low buffer consumption
- Multiple tray sizes and combs available

| Cat # | Tank Dimensions [W x L x H] | Gel Dimensions [mm] | Max. Sample | Capacity [mm] | Rapid Casting Gel |
|---------|----------------------------------|---------------------|-----------------|-------------------|-------------------|
| 03-3070 | 210 x 90 x 90mm | 70 x 70, 70 x 100 | 32 for 70 x 70 | 64 for 70 x 100 | Offer gel maker |
| 03-3100 | 220 x 125 x 90mm | 100 x 70, 100 x 100 | 50 for 100 x 70 | 100 for 100 x 100 | Offer gel maker |

BluView Transilluminator

Super bright blue LED lights can be used with a variety of safe dyes, and the lights are harmless and safe to view the experiment result. You can cut the gel with accuracy. With its compact size and light weight, the transilluminator occupies little of your lab bench and is easy to move from lab to lab.



- Compact, and easy to carry
- Aluminum alloy casing design
- Low heat dispersion
- Energy saving product
- 470nm harmless blue light for direct human contact
- Thinner and more lightweight body
- Use accessories to document experiment results by a smart phone

| Cat # | Blue Light Wavelength [nm] | Base Design | Construction Material | Power Adaptor | Filter Viewing Surface Dimensions [mm] |
|---------|----------------------------|-------------|-----------------------|--------------------|---|
| 03-6210 | 470 | Flat bed | Aluminum alloy | 12VDC/5A, 60W Max. | 153×153 |


The GelSMART gel imaging system

Features

- Compact foot print with a field of 10x10cm
- Simple operation of experiments and quick acquisition of results
- Epi-blue light diminishes UV exposure to lab personnel and DNA damage in the sample.
- Adjustable epi-blue light source emits light with high intensity and minimum heat for better light control.
- Filter holders can accommodate most popular smartphones/tablets smartphones/tablets.
- Optional light-weight app



Applications



| Cat # | Туре | Light Source | White Light Plate | Filter Protective Shield | Power Supply |
|---------|---------------------------------|-------------------------|----------------------|---|--------------------------------|
| 01-420X | The GelSMART gel imaging system | 470nm epi blue light | 10-400/1 | Amber, using when observing or cutting gels | 12V / 1.5A AC power adapter |

Application Guidelines

Maximum centrifugal acceleration force can be influenced by following criteria:

- Use of organic solutions, temperature below 0°C, centrifuging without centrifuge adapter, etc. Conduct a trial run to determine proper conditions before any critical use begins.
- The RCF (relative centrifugal force) or G ratings for our disposable centrifuge tubes have been established at room temperature using tubes filled to nominal capacity with water and centrifuging for 30 minutes. recommended carriers, adapters, and cushions that fully support the tubes. If an angle head rotor is used or proper support is not provided, G values will be lower.
- Do not confuse speed or revolutions per minute (RPM) with relative centrifugal force (RCF). Instructions for centrifuging a sample at a given RPM and time are incomplete unless the rotor or radius is specified. Protocols should always state the time and RCF value for centrifuging a sample.

Autoclaving guideline for Plastics

To ensure successful steam sterilization of your plastic labware please adhere to the following guidelines:

PLEASE NOTE THAT NOT ALL PRODUCTS ARE AUTOCLAVABLE. Autoclave plastic consumables at 121°C and 15psi/1 atm for 15 minutes. Do not use a dry cycle, as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to 80°C or below before removing products.

Our Commitment to Quality

At Biologix, we look at quality from every angle in terms of functionality, safety, purity, formulation, consistency, manufacturing, customer service, technical support, and much more, all in order to consistently ensure customer satisfaction.

According to the detailed QA/QC procedures for our manufacturing facility, each and every product is traceable by lot, regardless of packaging format. Products must go through a five-step human inspection before being released for packaging.

Our facility is equipped with the latest instruments and technology from the most trusted manufacturers, ensuring that Biologix applies the most stringent standards to certify that products are free from RNase, DNase, endotoxin (pyrogens), etc.

In addition to being ISO 9001:2015, CE, and FDA compliant, we also make available to our customers lot-specific certificates such as Sterilization Certificates, Nuclease Free Certificates, Certificates of Origin, etc. Please contact Biologix for Certificates of Compliance.

Note: Product parameters within this catalogue are for reference only.

Chemical Resistance Form

For Biologix Plastic Laboratory Supplies

| | PS | PP | HDPE | LDPE | PVC | CA | PC | CN | NY | MCE | PTFE | PET |
|----------------------------------|----|----|------|------|-----|----|----|----|----|-----|------|-----|
| Acids | | | | | | | | | | | | |
| Hydrochloric acid (25%) | G | G | G | G | G | N | R | R | N | 0 | R | R |
| Hydrochloric acid (concentrated) | F | G | G | G | F | N | R | N | N | N | R | 0 |
| Nitric acid (concentrated) | Р | Р | N | N | Р | N | R | N | N | N | 0 | N |
| Nitric acid (25%) | Р | G | F | F | F | N | R | L | N | 0 | R | R |
| Alcohols | | | | | | | | | | | | |
| Butanol | G | G | G | G | G | R | R | R | R | R | R | R |
| Ethanol | G | G | G | G | G | R | R | N | R | 0 | R | R |
| Methanol | G | G | G | G | G | R | R | N | R | 0 | R | R |
| Amines | | | | | | | | | | | | |
| Aniline | G | G | G | G | Р | N | N | R | R | N | R | 0 |
| Dimethylformamide | Р | G | G | G | F | N | N | N | R | N | R | N |
| Bases | | | | | | | | | | | | |
| Ammonium hydroxide (25%) | F | G | G | G | G | R | N | R | R | 0 | N | 0 |
| Ammonium hydroxide (1N) | F | G | G | G | G | N | N | R | R | 0 | N | N |
| Sodium hydroxide | G | G | G | G | G | N | N | N | R | N | R | N |
| Hydrocarbons | | | | | | | | | | | | |
| Hexane | Р | G | G | G | F | R | R | R | R | R | R | R |
| Toluene | Р | G | N | N | Р | R | 0 | R | R | R | R | N |
| Xylene | Р | F | Р | N | Р | R | R | R | R | R | R | N |
| Dioxane | Р | G | G | F | Р | N | N | N | R | N | R | R |
| Dimethylsulfoxide (DMSO) | Р | G | G | G | Р | N | N | N | R | N | R | 0* |
| Halogenated Hydrocarbons | | | | | | | | | | | | |
| Chloroform | Р | G | Р | N | Р | N | N | R | R | N | R | R |
| Methylene chloride | Р | F | Р | N | Р | N | N | R | R | N | R | N |
| Ketones | | | | | | | | | | | | |
| Acetone | Р | G | G | G | Р | N | 0 | N | R | N | R | R |
| Methyl ethyl diketone | Р | G | G | G | Р | N | 0 | N | R | 0 | R | R |

^{*} Can be used with aqueous solutions containing up to 20% DMSO.

The above chemical resistance chart is meant to be used as a guide only.

LEGEND

N = Not Recommended P = Poor PVC = Polyvinyl Chloride NY = Nylon

Biologix Group Limited

Phone: +86-531-67802668

Biologix Corporation

No.2766 Ying Xiu Road High-Tech Industrial Development Zone Jinan, Shandong 250101 P.R. China

Biologix Plastics (Changzhou) Co.,Ltd.

No.45 Chang Jin East Road, Zouqu Town, Zhonglou District, Changzhou City, Jiangsu, 213141, P. R. China

Biologix USA, Inc.

931 Via Alondra, Camarillo, CA 93012, United States

Phone: 1-805 586 6718

Email: Customercare@Biologixusa.com

www.BiologixUSA.com

Biologix Europe GmbH

Ludwigstrasse 44, 85399 Hallbergmoos, Germany

Email: Order@Biologixeurope.com

www.BiologixEurope.com





Inquiry@BiologixGroup.com www.BiologixGroup.com www.CryoKING.com